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
WESTERN BLOT APPLICATION GUIDE

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Expressionⁱ

Tissue specificityⁱ

Expressed in prostate cancer and levels increase from the normal to the malignant state (at protein level). Expressed in all human cell types so far analyzed. The Tyr-176 phosphorylated form shows a significant increase in expression in breast cancers during the progressive stages i.e. normal to hyperplasia (ADH), ductal carcinoma in situ (DCIS), invasive ductal carcinoma (IDC) and lymph node metastatic (LNMM) stages. [3 Publications](#)

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PTM/Processingⁱ

Features

Showing features for chain¹, modified residue¹, modified residue (large scale data)¹, disulfide bond¹, glycosylation¹, cross-link¹.

[Download](#)

TYPE	ID	POSITION(S)	SOURCE	DESCRIPTION
Modified residue		126	UniProt	Phosphoserine; alternate Combined Sources
Glycosylation		126	UniProt	O-linked (GlcNAc) serine; alternate 1 Publication
Modified residue (large scale data)		126	PRIDE	Phosphoserine Combined Sources
Modified residue		129	UniProt	Phosphoserine; alternate Combined Sources
Glycosylation		129	UniProt	O-linked (GlcNAc) serine; alternate 1 Publication
Modified residue (large scale data)		129	PRIDE	Phosphoserine Combined Sources
Modified residue (large scale data)		140	PTMeXchange	Ubiquitinated lysine Silver Combined Sources
Modified residue (large scale data)		154	PTMeXchange	Ubiquitinated lysine Gold Combined Sources
Modified residue (large scale data)		158	PTMeXchange	Ubiquitinated lysine Silver Combined Sources

[Expand Table](#)

Post-translational modificationⁱ

O-GlcNAcylation at Thr-305 and Thr-312 inhibits activating phosphorylation at Thr-308 via disrupting the interaction between AKT1 and PDPK1. O-GlcNAcylation at Ser-473 also probably interferes with phosphorylation at this site. [14 Publications](#)

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Align isoforms (2) Add isoforms

Sequence status¹ Complete

This entry describes 2 isoforms¹ produced by Alternative splicing.

P31749-1

This isoform has been chosen as the canonical sequence. All positional information in this entry refers to it. This is also the sequence that appears in the downloadable versions of the entry.

Name 1

See also sequence in UniParc or sequence clusters in UniRef

Tools Download Add Highlight Copy sequence

Length 400
Mass (Da) 55,686

Last updated 2005-02-01 v2
MD5 Checksum¹ 1620296E269940819BAFA78CA71E003A

MSDVAIVEKGLMLHKRGEYIKTNRPRYFLKNDGTFIGYKERPQDVQREAPLNNF5VAQCQLNKTERTPRNTFIIRCLQNTTVIERTFHWETPEEREENTTAIQTVDGLKKQEEEEPDFRSGSPSDNSGAEEHEVSL
PKHRVTNMEFEYLKLLGKGTFGKVLVKEKATGRYAMKILKKEVIVAKDEVHTLTENRVLQNSRHPFLTALKYSFQTHDRLCFVHEYANGGELFFHLSRERVFSEDRARFYGAEIVSALDYLHSEKMYVYRDLKLE
MLDKDGHGIIKIDTDFGLCKEGIKDGMTHKTFCTGPEYLAPEVEDNDYGRAVDMMGLGVVHYEMMKGRLPFYNQDHEKLFELILMEETRFPRTLGPEAKSLLSGLLKDKPKQRLGGGSEDAKEIMQHRFFAGIVWQHVE
LSPPFKPQVTSETDTRYFDEEFTAQHITITPPDQDSMECVDSERRPHFPQFSYASAGTA

P31749-2

Name 2

Differences from canonical 1-62: Missing 1 Publication

See also sequence in UniParc or sequence clusters in UniRef

Show sequence

Computationally mapped potential isoform sequences¹

There are 6 potential isoforms mapped to this entry

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
RNA EXPRESSION OVERVIEW

Organ Origin Category Expression Alphabetical

Cell lines ordered by descending RNA expression order.

Cell Line	Organ	Origin	Category	Expression (nTPM)
ASC2xelo	Brain			23
Sfha	Brain			20
K-562	Liver & Gallbladder			19
A-549	Gastrointestinal tract			18
HaCat	Male reproductive system			17
HTERT-RE1	Kidney & Urinary bladder			16
Rh30	Skin			15
hep-G2	Eye			14
A272	Proximal digestive tract			13
OE19	Lung			12
JURKAT	Female reproductive system			11
HOLM-2	Endothelial			10
PC-3	Mesenchymal			9
REL	Lymphoid			8
CAO-2	Myeloid			7
HTCEP	Renal, urinary, male reproductive system			6
NBA	Brain			5
B/J Human fibroblast	Brain			4
PDON TERT26	Brain			3
HBA293	Brain			2
A-431	Brain			1
GMG	Brain			1
MCF-1	Brain			1
U-251MG	Brain			1
HTERT-HNE1	Brain			1
SUSa	Brain			1
HeLa	Brain			1
REL	Brain			1
RPTCE TERT1	Brain			1
HTEC-15TER 24.8	Brain			1
HBB TERT88	Brain			1
U2OS	Brain			1
EFO21	Brain			1
Thp-1	Brain			1
HDFTERT166	Brain			1
SK-MEL-30	Brain			1
HAP1	Brain			1
SK-SY5Y	Brain			1
HBCRC-MT	Brain			1
PDON SW TERT153	Brain			1
HUVEC TERT2	Brain			1
IF-4	Brain			1
LYCN-MO2	Brain			1

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WESTERN BLOT APPLICATION GUIDE

产品类别	货号	产品名称
预制胶、快速制胶试剂盒、或者散装配胶等	BRK0003	一步法彩色快速凝胶试剂盒 10%
	BRK0004	一步法彩色快速凝胶试剂盒 12%
	BRK0005	一步法彩色快速凝胶试剂盒 15%
	BRK1010T	预制 Tris-Gly (10%) 10 孔超快蛋白电泳试剂盒
	BRK1210T	预制 Tris-Gly (12%) 10 孔超快蛋白电泳试剂盒
	BRK0021P	预制 Tris-Gly (4-20%) 10 孔超快蛋白电泳试剂盒
	BRK1010B	预制 Bis-Tris (10%) 10 孔蛋白电泳试剂盒
	BRK1210B	预制 Bis-Tris (12%) 10 孔蛋白电泳试剂盒
	BRK0031	预制 Bis-Tris (4-20%) 10 孔蛋白电泳试剂盒
	BRK0020	预制 Tris-Gly (4-20%) 蛋白电泳试剂盒
	BRK0020P	预制 Tris-Gly (4-20%) 超快蛋白电泳试剂盒
	BRK0030	预制 Bis-Tris (4-20%) 蛋白电泳试剂盒
	BR00013	Cassette Opener (预制胶开胶器)
电泳、转膜、剥离液等试剂及耗材	RM00010	1X Tris-Glycine SDS Running Buffer
	BR00007	超快电泳缓冲液 (Tris- 甘氨酸 SDS 粉末)
	BR00009	Ice-Bath Free Ultra-Fast Transfer Buffer Powder { 免冰浴超快转膜缓冲液粉末 }
	BR00016	WB Stripping buffer 抗体洗脱液 (再生液、剥离液)
	BR00017	ecolorization-free Ultra Fast Staining Buffer 考马斯亮蓝快速染色液 (免脱色)
	RM00017	Nitrocellulose membrane (6.3×8.3cm)(0.45μm) {NC 膜 (0.45μm)}
	RM02801	Nitrocellulose membrane (6.3×8.3cm)(0.2μm) {NC 膜 (0.2μm)}
	RM00018	Polyvinylidene fluoride (6.3×8.3cm)(0.2μm) {PVDF 膜 }
洗膜、封闭、一抗二抗孵育及曝光	RM02954	Western Blot Primary Antibody Dilution Buffer { WB 一抗稀释液 }
	RM02955	Universal Antibody Dilution Buffer { 通用型抗体稀释液 }
	RM02956	Rapid Blocking Buffer { 快速封闭液 (5*) }
	RM00020P	ECL Basic Plus Kit {ECL 化学发光底物 Plus(普通型)}
	RM00021P	ECL Enhanced Plus Kit {ECL 化学发光底物 Plus(增强型)}
	RM02867	ECL Super Kit {ECL 化学发光底物 (超敏型)}
常见内参抗体	AC012	α-Tubulin Mouse mAb
	AC026	β-Actin Rabbit mAb (High Dilution)
	AC021	β-Tubulin Mouse mAb
	AC038	β-Actin Rabbit mAb
	A12289	β-Tubulin Rabbit mAb
	AC009	Actin (plant specific) Mouse mAb
	AC002	GAPDH Mouse mAb
	A19056	GAPDH Rabbit mAb(High Dilution)
	AC033	GAPDH Mouse mAb (High Dilution)
	AC028	HRP-conjugated β-Actin Rabbit mAb

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PMID	标题	期刊名称	影响因子	产品类型	引用货号/物种
25083871	MicroRNA directly enhances mitochondrial translation during muscle differentiation	Cell	41.582	抗体	A4967
28753426	Methyltransferase SETD2-Mediated Methylation of STAT1 Is Critical for Interferon Antiviral Activity	Cell	41.582	抗体	A3194
28965759	Temporal Control of Mammalian Cortical Neurogenesis by m6A Methylation	Cell	41.582	抗体	A8530
29290468	Stabilization of Reversed Replication Forks by Telomerase Drives Telomere Catastrophe	Cell	41.582	抗体	A9555
30146159	Microbial Siderophore Enterobactin Promotes Mitochondrial Iron Uptake and Development of the Host via Interaction with ATP Synthase.	Cell	41.582	抗体	A5769
31031002	Structure and Degradation of Circular RNAs Regulate PKR Activation in Innate Immunity.	Cell	41.582	抗体	AP0342
31251911	Pervasive Chromatin-RNA Binding Protein Interactions Enable RNA-Based Regulation of Transcription.	Cell	41.582	抗体	A5016
31585088	Integrated Proteogenomic Characterization of HBV-Related Hepatocellular Carcinoma	Cell	41.582	抗体	A11932
32169215	Initiation of Parental Genome Reprogramming in Fertilized Oocyte by Splicing Kinase SRPK1-Catalyzed Protamine Phosphorylation	Cell	41.582	抗体	A5854
32259487	Distinct Processing of lncRNAs Contributes to Non-conserved Functions in Stem Cells	Cell	41.582	抗体	A5964
32810437	Formation of NPR1 Condensates Promotes Cell Survival during the Plant Immune Response	Cell	41.582	抗体	A20262
32931733	Targeting Mitochondria-Located circRNA SCAR Alleviates NASH via Reducing mROS Output	Cell	41.582	抗体	A9840
33991488	Mouse totipotent stem cells captured and maintained through spliceosomal repression	Cell	41.58	抗体	A5891, A6709
28297716	RNA m6A methylation regulates the ultraviolet-induced DNA damage response.	Nature	49.962	抗体	A1848, AP0099
30111841	Modulating plant growth-metabolism coordination for sustainable agriculture	Nature	49.962	抗体	A16279
30814741	Proteomics identifies new therapeutic targets of early-stage hepatocellular carcinoma.	Nature	49.962	抗体	A2091
31485071	Structural basis of nucleosome recognition and modification by MLL methyltransferases.	Nature	49.962	抗体	A2355, A2357
32376953	Plant 22-nt siRNAs mediate translational repression and stress adaptation	Nature	49.962	抗体	A2348
32238924	U1 snRNP regulates chromatin retention of noncoding RNAs	Nature	49.962	抗体	A5761, AP0844
32581380	Reversing a Model of Parkinson's Disease With in Situ Converted Nigral Neurons	Nature	49.962	抗体	A14150, A6107
33268895	Galactosaminogalactan activates the inflammasome to provide host protection	Nature	49.962	抗体	A15094, A6724
33361816	Molecular basis of nucleosomal H3K36 methylation by NSD methyltransferases	Nature	49.962	抗体	A2364
32029600	Enhanced sustainable green revolution yield via nitrogen-responsive chromatin modulation in rice.	Science	47.728	抗体	A16279

WESTERN BLOT APPLICATION GUIDE

PMID	标题	期刊名称	影响因子	产品类型	引用货号/物种
33243860	Identification of Integrator-PP2A complex (INTAC), an RNA polymerase II phosphatase	Science	47.728	抗体	A2355, A7253
33795473	Structural insights into preinitiation complex assembly on core promoters	Science	47.728	抗体	A2192, A5928
34709882	NIN-like protein transcription factors regulate leghemoglobin genes in legume nodules	Science	47.72	抗体	A20191
30348677	Spatiotemporal loss of NF1 in Schwann cell lineage leads to different types of cutaneous neurofibroma susceptible to modification by the Hippo pathway.	Cancer discovery	39.397	抗体	A8202
33653692	Macropinocytosis in Cancer-Associated Fibroblasts is Dependent on CaMKK2/ARHGEF2 Signaling and Functions to Support Tumor and Stromal Cell Fitness	Cancer Discov	39.39	抗体	A16678
33558758	Comprehensive characterization of protein-protein interactions perturbed by disease mutations	Nat Genet	38.33	抗体	A2877
28263319	Asymmetric subgenome selection and cis-regulatory divergence during cotton domestication	Nature genetics	38.33	抗体	A2363
30193849	The Zinc-Finger Protein ZCCHC3 Binds RNA and Facilitates Viral RNA Sensing and Activation of the RIG-I-like Receptors.	Immunity	31.745	抗体	A12938
30650373	Fever Promotes T Lymphocyte Trafficking via a Thermal Sensory Pathway Involving Heat Shock Protein 90 and α 4 Integrins.	Immunity	31.745	抗体	A0284, A0969
30770248	Specific Decrease in B-Cell-Derived Extracellular Vesicles Enhances Post-Chemotherapeutic CD8+ T Cell Responses.	Immunity	31.745	抗体	A10389, A1934
29033244	iASPP Is an Antioxidative Factor and Drives Cancer Growth and Drug Resistance by Competing with Nrf2 for Keap1 Binding	Cancer Cell	31.743	抗体	A0207, A0208
32049046	Direct Phosphorylation and Stabilization of MYC by Aurora B Kinase Promote T-cell Leukemogenesis	Cancer Cell	31.743	抗体	A11021, A5904
33186520	Elevated CXorf67 Expression in PFA Ependymomas Suppresses DNA Repair and Sensitizes to PARP Inhibitors	Cancer Cell	31.743	抗体	A2189
32442403	SETD5-Coordinated Chromatin Reprogramming Regulates Adaptive Resistance to Targeted Pancreatic Cancer Therapy	Cancer Cell	31.743	抗体	A7254, A7257
33740421	Cancer cells escape autophagy inhibition via NRF2-induced macropinocytosis	Cancer Cell	31.74	抗体	A11159
31915373	Single-cell reconstruction of the adult human heart during heart failure and recovery reveals the cellular landscape underlying cardiac function	Nature Cell Biology	28.824	抗体	A2179
32029896	BCAT2-mediated BCAA catabolism is critical for development of pancreatic ductal adenocarcinoma	Nature Cell Biology	28.824	抗体	A1957
29510730	The lncRNA MACC1-AS1 promotes gastric cancer cell metabolic plasticity via AMPK/Lin28 mediated mRNA stability of MACC1	Molecular Cancer	27.401	抗体	A0207, A2156
29803224	The long non-coding RNA PTTG3P promotes cell growth and metastasis via up-regulating PTTG1 and activating PI3K/AKT signaling in hepatocellular carcinoma	Molecular Cancer	27.401	抗体	A0266
30458784	Invasion-related circular RNA circFND3B inhibits bladder cancer progression through the miR-1178-3p/G3BP2/SRC/FAK axis.	Molecular cancer	27.401	抗体	AP0185, AP0302
30541550	TRIB2 functions as novel oncogene in colorectal cancer by blocking cellular senescence through AP4/p21 signaling.	Molecular cancer	27.401	抗体	A11661
30922314	RNA N6-methyladenosine demethylase FTO promotes breast tumor progression through inhibiting BNIP3.	Molecular cancer	27.401	抗体	A0208, A2156
30940151	Circular RNA circTADA2A promotes osteosarcoma progression and metastasis by sponging miR-203a-3p and regulating CREB3 expression.	Molecular cancer	27.401	抗体	A0289
31122259	miR-301a promotes lung tumorigenesis by suppressing Runx3.	Molecular cancer	27.401	抗体	A0433, A6247
31718709	Circ-HuR suppresses HuR expression and gastric cancer progression by inhibiting CNBP transactivation	Molecular Cancer	27.401	抗体	A0958

WESTERN BLOT APPLICATION GUIDE

PMID	标题	期刊名称	影响因子	产品类型	引用货号/物种
29275960	Inactivation of the Glucose-Dependent Insulinotropic Polypeptide Receptor Improves Outcomes following Experimental Myocardial Infarction	Cell Metabolism	27.287	抗体	A9816
30527742	The Circadian Protein Period2 Suppresses mTORC1 Activity via Recruiting Tsc1 to mTORC1 Complex.	Cell Metabolism	27.287	抗体	A2190, A3217
31230984	OLFR734 Mediates Glucose Metabolism as a Receptor of Asprosin.	Cell Metabolism	27.287	抗体	A2036
31668872	Adipocyte Hypoxia-Inducible Factor 2 α Suppresses Atherosclerosis by Promoting Adipose Ceramide Catabolism	Cell Metabolism	27.287	抗体	A10613
31761563	Fructose-1,6-Bisphosphatase 2 Inhibits Sarcoma Progression by Restraining Mitochondrial Biogenesis	Cell Metabolism	27.287	抗体	A13552, A6774
32375062	Low-Dose Sorafenib Acts as a Mitochondrial Uncoupler and Ameliorates Nonalcoholic Steatohepatitis	Cell Metabolism	27.287	抗体	A11435, AP0115
33186558	Elevation of JAML Promotes Diabetic Kidney Disease by Modulating Podocyte Lipid Metabolism	Cell Metabolism	27.287	抗体	A13049
30174306	Short-Term Mitochondrial Permeability Transition Pore Opening Modulates Histone Lysine Methylation at the Early Phase of Somatic Cell Reprogramming	Cell Metabolism	27.287	抗体	A1156, A2359
30455424	VCAM-1+ macrophages guide the homing of HSPCs to a vascular niche.	Nature	49.962	抗体定制	斑马鱼
30944466	TMK1-mediated auxin signalling regulates differential growth of the apical hook.	Nature	49.962	抗体定制	拟南芥
31942069	Host-mediated ubiquitination of a mycobacterial protein suppresses immunity	Nature	49.962	抗体定制	结核杆菌
30487604	Stella safeguards the oocyte methylome by preventing de novo methylation mediated by DNMT1.	Nature	49.962	抗体定制	小鼠
30356214	Nuclear cGAS suppresses DNA repair and promotes tumorigenesis.	Nature	49.962	抗体定制	人
31243371	UDP-glucose accelerates SNAIL mRNA decay and impairs lung cancer metastasis.	Nature	49.962	抗体定制	人
32841601	A Defense Pathway Linking Plasma Membrane and Chloroplasts and Co-opted by Pathogens	Cell	41.582	抗体定制	拟南芥
28552346	Ubiquitination-Deficient Mutations in Human Piwi Cause Male Infertility by Impairing Histone-to-Protamine Exchange during Spermiogenesis	Cell	41.582	抗体定制	小鼠
31835033	A Translation-Activating Function of MIWI/piRNA during Mouse Spermiogenesis	Cell	41.582	抗体定制	小鼠
30612740	METTL13 Methylation of eEF1A Increases Translational Output to Promote Tumorigenesis.	Cell	41.582	抗体定制	人
31859049	Immune Sensing of Cell Death through Recognition of Histone Sequences by C-Type Lectin-Receptor-2d Causes Inflammation and Tissue Injury	Immunity	31.745	抗体定制	小鼠
31085179	Inhibition of Nuclear PTEN Tyrosine Phosphorylation Enhances Glioma Radiation Sensitivity through Attenuated DNA Repair	Cancer Cell	31.743	抗体定制	人和小鼠
27002220	A positive role for polycomb in transcriptional regulation via H4K20me1.	Cell research	25.617	抗体定制	果蝇
30778176	An H3K27me3 demethylase-HSFA2 regulatory loop orchestrates transgenerational thermomemory in Arabidopsis.	Cell research	25.617	抗体定制	拟南芥
25753158	ERK kinase phosphorylates and destabilizes the tumor suppressor FBW7 in pancreatic cancer.	Cell research	25.617	抗体定制	人
26928339	The kinase CK1 ϵ controls the antiviral immune response by phosphorylating the signaling adaptor TRAF3.	Nature immunology	25.606	抗体定制	小鼠
32461121	Peptide SMIM30 promotes HCC development by inducing SRC/YES1 membrane anchoring and MAPK pathway activation	J Hepatol	25.083	抗体定制	小鼠

WESTERN BLOT APPLICATION GUIDE

PMID	标题	期刊名称	影响因子	产品类型	引用货号/物种
32735848	A Bacterial Effector Protein Hijacks Plant Metabolism to Support Pathogen Nutrition	Cell Host Microbe	21.023	抗体定制	烟草
32994405	The cross-talk between methylation and phosphorylation in lymphoid-specific helicase drives cancer stem-like properties	Signal Transduct Target Ther	18.187	抗体定制	人
32296027	CDK4/6 inhibition blocks cancer metastasis through a USP51-ZEB1-dependent deubiquitination mechanism	Signal Transduction and Targeted Therapy	18.187	抗体定制	人
31422874	Naa10p Inhibits Beige Adipocyte-Mediated Thermogenesis through N- α -acetylation of Pgc1 α .	Molecular cell	17.97	抗体定制	小鼠
31110366	Oxidization of TGF β -activated kinase by MPT53 is required for immunity to Mycobacterium tuberculosis.	Nature microbiology	17.745	抗体定制	小鼠
33619562	The human tRNA taurine modification enzyme GTPBP3 is an active GTPase linked to mitochondrial diseases.	Nucleic Acids Res	16.971	抗体定制	人
33176558	A VAST-domain protein regulates autophagy, membrane tension, and sterol homeostasis in rice blast fungus	Autophagy	16.016	抗体定制	稻瘟病菌
25814387	Suppressor of Deltex mediates Pez degradation and modulates Drosophila midgut homeostasis.	Nat Commun	14.919	抗体定制	果蝇
30833710	In rose, transcription factor PTM balances growth and drought survival via PIP2;1 aquaporin	Nature plant	13.256	抗体定制	月季
31075443	Translational Regulation of Plant Response to High Temperature by a Dual-Function tRNA His Guanylyltransferase in Rice	Molecular Plant	13.164	抗体定制	水稻
30914468	OsSHI1 Regulates Plant Architecture Through Modulating the Transcriptional Activity of IPA1 in Rice	Plant cell	11.277	抗体定制	水稻
30923230	Rice qGL3/OsPPKL1 Functions with the GSK3/SHAGGY-Like Kinase OsGSK3 to Modulate Brassinosteroid Signaling	Plant cell	11.277	抗体定制	水稻
31530735	Intra-Kernel Reallocation of Proteins in Maize Depends on VP1-Mediated Scutellum Development and Nutrient Assimilation	Plant cell	11.277	抗体定制	玉米
30962394	Maize VKS1 Regulates Mitosis and Cytokinesis during Early Endosperm Development	Plant cell	11.277	抗体定制	玉米
32155285	miR-210 participates in hepatic ischemia reperfusion injury by forming a negative feedback loop with SMAD4	Hepatology (Baltimore, Md.)	17.425	ELISA	RK00006
33305325	Entropy subspace separation-based clustering for noise reduction (ENCORE) of scRNA-seq data	Nucleic Acids Res	16.971	ELISA	RK00719
33009518	Targeting a positive regulatory loop in the tumor-macrophage interaction impairs the progression of clear cell renal cell carcinoma	Cell Death Differ	15.828	ELISA	RK00077
31924786	Arf1-mediated lipid metabolism sustains cancer cells and its ablation induces anti-tumor immune responses in mice	Nature Communications	14.919	ELISA	RK01574
34678143	Glycogen accumulation and phase separation drives liver tumor initiation	Cell	41.582	NGS	RK20351,RK20303
33523829	Abnormal neocortex arealization and Sotos-like syndrome-associated behavior in Setd2 mutant mice	Sci Adv	14.136	NGS	RK20202,RK20282
31585088	Integrated Proteogenomic Characterization of HBV-Related Hepatocellular Carcinoma	Cell	41.582	分子酶	RK20400
32017896	LARP7-Mediated U6 snRNA Modification Ensures Splicing Fidelity and Spermatogenesis in Mice	Molecular Cell	17.97	分子酶	RK20400
33473107	Endonuclease G promotes autophagy by suppressing mTOR signaling and activating the DNA damage response	Nat Commun	14.919	分子酶	RK20400,RK21203
32867860	Bradymonabacteria, a novel bacterial predator group with versatile survival strategies in saline environments	Microbiome	14.65	分子酶	RK21203
32789178	NuRD mediates mitochondrial stress-induced longevity via chromatin remodeling in response to acetyl-CoA level	Sci Adv	14.136	分子酶	RK21020
33523838	N6-methyldeoxyadenine and histone methylation mediate transgenerational survival advantages induced by hormetic heat stress	Sci Adv	14.136	分子酶	RK20400,RK21203

为全球生命科学工作者提供更好的产品和服务
Empowering Global Life Science with Superior Products and Services

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