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## Cell lysis buffer for IP (without inhibitors)

**Catalog number:** RM00022

**Volume:** 100 mL

### Product Introduction

This Buffer is a mammalian whole cell lysis buffer based on a modified RIPA buffer formulation (the Ingredients are shown as follows). This Lysis Buffer is compatible with BCA Protein Assay or Bradford Protein Assay. The extracted protein samples with this lysis buffer can be used in Western Blot, IP and Co-IP, etc. This lysis Buffer does not contain protease or phosphatase inhibitors. If desired, add Protease and Phosphatase Inhibitor Cocktail immediately before use.

### Reagents of Buffer

(Ingredients and final Concentration)

NaCl 136.89mM; Na<sub>2</sub>HPO<sub>4</sub> 1.01mM; KH<sub>2</sub>PO<sub>4</sub> 1.76mM; KCl 2.68 mM; SDS 0.347 mM

NP-40 0.5% (V/V)

### Storage

Can be stored at - 20°C within 12 months.

### Product Usage Information

Thaw this lysis buffer at 24-30°C, mixing end-over-end.

Aliquot proper amount of buffer based on cell numbers and add Protease and Phosphatase Inhibitor Cocktail immediately before use.

For lysis of adherent cells, we recommend the following: (all reagents and lysates must be kept cold)

1. Carefully wash cultured cells with pre-chilled PBS for 2 times.
2. Add in pre-chilled cell lysis buffer with Protease and Phosphatase Inhibitor Cocktail (1mL for 10<sup>7</sup> cells, 100~200μL for one well of 6-well plate).
3. Scrap cells off to clean 1.5mL tubes with a clean, cold scraper. Put them on a low-speed rotating shaker for 15 min at 4°C.
4. Centrifuge at 12,000 g 4°C for 5min, transfer the supernatant to new tubes immediately.
5. The supernatants are extracted protein samples, which can be used in Western Blot, IP and Co-IP, etc. after measuring protein concentration.

For lysis of non-adherent cells, we recommend the following: (all reagents and lysates must be kept cold)

1. Collect cell pellets by centrifugation (1000g) and remove residual culture media.
2. Carefully wash cells by resuspending cells with pre-chilled PBS.
3. Collect cell pellets by centrifugation (1000g) and remove residual PBS.
4. Repeat step 2 and 3.
5. Add in pre-chilled cell lysis buffer pre-chilled cell lysis buffer with Protease and Phosphatase Inhibitor Cocktail (1mL for 10<sup>7</sup> cells) and mix by pipetting.
6. Put them on a low-speed rotating shaker for 15 min at 4°C.
7. Centrifuge at 12,000 g 4°C for 5min, transfer the supernatant to new tubes immediately.
8. The supernatants are extracted protein samples, which can be used in Western Blot, IP and Co-IP, etc. after measuring protein concentration.

For lysis of tissue samples, we recommend the following: (all reagents and lysates must be kept cold)

1. Place tissue sample in a 60 mm or 100 mm dish and finely mince using a clean scalpel or razor blade. Keep dish on ice. It is important to keep the tissue cold to avoid protein degradation.
2. Collect the tissue sample pieces in 1.5mL tube and add pre-chilled cell lysis buffer with Protease and Phosphatase Inhibitor Cocktail (100~200 $\mu$ L for 20mg tissue).
3. Homogenize the tissue lysates by Dounce Homogenizer or other Homogenizers on ice as soon as possible.
4. Centrifuge at 12,000 g 4°C for 5min, transfer the supernatant to new tubes immediately. The supernatants are extracted protein samples, which can be used in Western Blot, IP and Co-IP, etc. after measuring protein concentration.