


Validation Data



Flow cytometry: 1×10^6 U-118 MG cells (negative control, left) and MCF7 cells (right) were surface-stained with ABflo® 594 Rabbit anti-Human TROP-2 mAb (A24177, 5 μ l/Test, orange line) or ABflo® 594 Rabbit IgG isotype control (A23821, 5 μ l/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).

Flow cytometry: 1×10^6 MCF7 cells were surface-stained with ABflo® 594 Rabbit IgG isotype control (A23821, 5 μ l/Test, left) or ABflo® 594 Rabbit anti-Human TROP-2 mAb (A24177, 5 μ l/Test, right).