


Validation Data



Flow cytometry: 1×10^6 HeLa cells (negative control, left) and SH-SY5Y cells (right) were surface-stained with Mouse anti-Human CD56 mAb (orange line) or PerCP Mouse IgG 2a isotype control (A26259,5 $\mu\text{l}/\text{Test}$, blue line). Non-fluorescently stained cells were used as blank control (red line).

Flow cytometry: 1×10^6 SH-SY5Y cells were surface-stained with PerCP Mouse IgG 2a isotype control (A26259,5 $\mu\text{l}/\text{Test}$, left) or Mouse anti-Human CD56 mAb (right).