



Validation Data



Flow cytometry: 1X10⁶ HeLa cells (negative control, left) and SH-SY5Y cells (right) were surface-stained with Mouse anti-Human CD56 mAb (orange line) or PerCP Mouse IgG 2a isotype control (A26259,5 µl/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).

Flow cytometry: 1X10⁶ SH-SY5Y cells were surface-stained with PerCP Mouse IgG 2a isotype control (A26259,5 µl/Test, left) or Mouse anti-Human CD56 mAb (right).