



Validation Data



Flow cytometry: 1X10⁶ NIH/3T3 cells (Low Expression, left) and RAW 264.7 cells (right) were intracellularly-stained with ABflo® 700 Rabbit anti-Mouse CD68 mAb (A26452, 5 µl/Test, orange line) or ABflo® 700 Rabbit IgG isotype control (A25976, 5 µl/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).

Flow cytometry: 1X10⁶ RAW 264.7 cells were intracellularly-stained with ABflo® 700 Rabbit IgG isotype control (A25976, 5 µl/Test, left) or ABflo® 700 Rabbit anti-Mouse CD68 mAb (A26452, 5 µl/Test, right).