


Validation Data



Flow cytometry: 1×10^6 NTERA-2 cells were intracellularly-stained with PE Rabbit IgG isotype control (A24172, 5 μ l/Test, left) or PE Rabbit anti-Human/Mouse Nanog mAb (A26535, 5 μ l/Test, right).

Flow cytometry: 1×10^6 F9 cells were intracellularly-stained with PE Rabbit IgG isotype control (A24172, 5 μ l/Test, left) or PE Rabbit anti-Human/Mouse Nanog mAb (A26535, 5 μ l/Test, right).