





## Validation Data



Flow cytometry: 1X10<sup>6</sup> HAP1 cells (Low Expression, left) and Daudi cells (right) were intracellularly-stained with ABflo® 610 Rabbit anti-Human/Mouse BCL6 mAb (A26765,5 µl/Test, orange line) or ABflo® 610 Rabbit IgG isotype control (A25826,5 µl/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).

Flow cytometry: 1X10<sup>6</sup> Daudi cells were intracellularly-stained with ABflo® 610 Rabbit IgG isotype control (A25826,5 µl/Test, left) or ABflo® 610 Rabbit anti-Human/Mouse BCL6 mAb (A26765,5 µl/Test, right).