


Validation Data



Flow cytometry: 1×10^6 CHO cells (negative control) and CHO(Transfection, right) cells were intracellularly-stained with PE Mouse anti DDDDK-Tag mAb(AE116, 5 μ l/Test, orange line) or PE Mouse IgG isotype control (5 μ l/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).

Flow cytometry: 1×10^6 CHO(Transfection) cells were intracellularly-stained with PE Mouse IgG isotype control (5 μ l/Test, left) or PE Mouse anti DDDDK-Tag mAb(AE116, 5 μ l/Test, right).