

Recombinant Human Tet2 Protein

Catalog No.: RP01590LQ **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	54790	Q6N021-1

Tags

N-Strep

Synonyms

MDS; KIAA1546;TET2;KIAA1546

Product Information

Source	Purification
HEK293 cells	> 90% by SDS-PAGE.

Endotoxin

Please contact us for more information.

Formulation

Supplied as a 0.22 µm filtered solution in
50mM HEPES,300mM NaCl,1mM
TCEP,10% Glycerol,pH 8.0.

Reconstitution

Background

Basic Information

Description

Recombinant Human Tet2 Protein is produced by HEK293 cells expression system.

Bio-Activity

Storage

This product is stable at $\leq -70^{\circ}\text{C}$ for up to 1 year from the date of receipt.
For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature.
Avoid repeated freeze/thaw cycles.

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn