



扫码立即参与活动

甄选千种兔单抗，订购即获**1.6倍**兑换膨胀金

线上商城优惠同享，高至**8倍积分** ✨ ✨

活动时间：2021年9月1日---2021年10月31日



满额有礼/惊喜享不停

货号	名称	应用	物种
A19654	[KO Validated] active + pro Caspase-3 Rabbit mAb	WB,IF	H,M,R
A17911	[KO Validated] APP Rabbit mAb	WB,IHC,IF,IP	H,M,R
A19677	[KO Validated] ATG5 Rabbit mAb	WB,IHC	H,M,R
A18642	[KO Validated] Bax Rabbit mAb	WB,IHC,IP	H,M
A19684	[KO Validated] Bax Rabbit mAb	WB,IHC,IF	H,M,R
A20227	[KO Validated] Bax Rabbit mAb	WB,IHC,IF,IP	H,M,R
A18676	[KO Validated] Caspase-9 Rabbit mAb	WB,IF	H,M,R
A11666	[KO Validated] CD168/RHAMM Rabbit mAb	WB,IHC	H,M
A19020	[KO Validated] CD44 Rabbit mAb	WB,IHC	H
A11136	[KO Validated] CDK4 Rabbit mAb	WB,IF	H
A19094	[KO Validated] CDKN1A/p21CIP1 Rabbit mAb	WB,IHC	H
A19095	[KO Validated] CDKN1B/p27KIP1 Rabbit mAb	WB,IHC	H,M,R
A11651	[KO Validated] CDKN2A/p16INK4a Rabbit mAb	WB,IHC,IF	H
A5025	[KO Validated] CDKN2A/p16INK4a Rabbit mAb	WB,IF	H
A19032	[KO Validated] c-Myc Rabbit mAb	WB,IF,ChIP	H,M
A19588	[KO Validated] CTCF Rabbit mAb	WB,IF,ChIP	H,M,R
A19037	[KO Validated] Cyclin B1 Rabbit mAb	WB,IHC,IP	H,M
A19038	[KO Validated] Cyclin D1 Rabbit mAb	WB,IHC,IF	H,M,R
A4912	[KO Validated] Cytochrome C Rabbit mAb	WB,IHC	H,M,R
A19679	[KO Validated] DNMT1 Rabbit mAb	WB,IF	H,M,R
A19577	[KO Validated] EZH2/KMT6 Rabbit mAb	WB,IHC	H,M
A19582	[KO Validated] Fas Rabbit mAb	WB,IHC	H,M,R
A11488	[KO Validated] FGF2 Rabbit mAb	WB,IHC,IF	H,M,R
A4443	[KO Validated] GRK2 Rabbit mAb	WB,IF	H,M,R
A19571	[KO Validated] HDAC1 Rabbit mAb	WB,IHC,IF	H,M,R
A0533	[KO Validated] HK1 Rabbit mAb	WB,IHC,IF,IP	H,M,R
A19529	[KO Validated] HMGB1 Rabbit mAb	WB,IHC,IF	H,M,R
A19622	[KO Validated] HuR/ELAVL1 Rabbit mAb	WB,IHC,IF,IP	H,M,R

货号	名称	应用	物种
A19524	[KO Validated] Lamin A/C Rabbit mAb	WB,IHC,IF,IP	H,M,R
A19665	[KO Validated] LC3B Rabbit mAb	WB,IHC,IF	H,M,R
A19079	[KO Validated] METTL3 Rabbit mAb	WB,IF	H,M
A19685	[KO Validated] MLKL Rabbit mAb	WB,IHC	H,M
A4624	[KO Validated] MTA2 Rabbit mAb	WB,IHC,IF	H
A19083	[KO Validated] N-Cadherin Rabbit mAb	WB,IHC,IF	H,M,R
A19653	[KO Validated] NF- κ B p65/RelA Rabbit mAb	WB,IHC,IP	H,M,R
A19596	[KO Validated] PARP1 Rabbit mAb	WB,IHC	H,M,R
A19135	[KO Validated] PD-L1/CD274 Rabbit mAb	WB,IHC	H,R
A19104	[KO Validated] PTEN Rabbit mAb	WB,IHC,IF	H,M,R
A3645	[KO Validated] RBBP4 Rabbit mAb	WB,IHC,IF	H,M,R
A19667	[KO Validated] SIRT1 Rabbit mAb	WB,IF	H
A19113	[KO Validated] Smad1 Rabbit mAb	WB,IF	H,M
A19114	[KO Validated] Smad2 Rabbit mAb	WB,IF	H,M,R
A19115	[KO Validated] Smad3 Rabbit mAb	WB,IHC,IF	H,M,R
A19116	[KO Validated] Smad4 Rabbit mAb	WB,IHC,IP	H,M,R
A3247	[KO Validated] SMARCB1/SNF5 Rabbit mAb	WB,IF	H,M,R
A19700	[KO Validated] SQSTM1/p62 Rabbit mAb	WB,IHC,IF,IP	H,M,R
A19119	[KO Validated] Src Rabbit mAb	WB,IF	H,M,R
A19566	[KO Validated] STAT3 Rabbit mAb	WB,IHC,IF,IP	H,M,R
A19666	[KO Validated] TTC11/FIS1 Rabbit mAb	WB,IHC,IF,IP	H,M,R
A19607	[KO Validated] Vimentin Rabbit mAb	WB,IHC,IF,IP	H,M,R
A19657	[KO Validated] β -Catenin Rabbit mAb	WB,IHC,IF	H,M,R
A9152	14-3-3 beta/zeta Rabbit mAb	WB,IF	H,M,R
A9162	14-3-3 gamma Rabbit mAb	WB,IHC	H,M,R
A4377	14-3-3 sigma Rabbit mAb	WB,IHC	H,M
A8936	14-3-3 Theta Rabbit mAb	WB,IF	H,M,R
A5024	15-PGDH/HPGD Rabbit mAb	WB,IF	M,R
A3859	53BP1 Rabbit mAb	WB,IHC,IF	H,M
A19706	5HT7 Receptor Rabbit mAb	WB,IF	H,M,R
A20599	5-Methylcytosine (5mC) Rabbit mAb	DB	H,M,R,Other (Wide Range)
A5968	67kDa Laminin Receptor Rabbit mAb	WB,IF	H,M,R
A9135	ABCE1 Rabbit mAb	WB,IF	H,M

货号	名称	应用	物种
A19621	ABCF1 Rabbit mAb	WB,IHC	H,R
A17907	ABCG1 Rabbit mAb	WB	H,M,R
A0664	ACAA2 Rabbit mAb	WB,IF	H,M,R
A4567	ACADM Rabbit mAb	WB,IHC,IF	H,M,R
A19627	ACC1 Rabbit mAb	WB,IHC	H,M,R
A11357	ACE1 Rabbit mAb	WB, IP	M,R
A4612	ACE2 Rabbit mAb	WB,IHC	H,M,R
A2770	Acetyl-Histone H3-K23 Rabbit mAb	WB,IF	H,M,R
A2771	Acetyl-Histone H3-K23 Rabbit mAb	WB,IHC,IF	H,M,R
A19525	Acetyl-Histone H4-K5 Rabbit mAb	WB,IHC,IF	H,M,R,Other (Wide Range)
A19836	Acetyl-p53- K370 Rabbit mAb	WB	H,M,R
A3719	ACLY Rabbit mAb	WB,IF	H,M,R
A4524	Aconitase 2 (ACO2) Rabbit mAb	WB,IHC,IF	H,M,R
A11021	Active Caspase-3 Rabbit mAb	WB	H
A7709	ACTL6A Rabbit mAb	WB,IHC	H,M,R
A8939	ACTN2 Rabbit mAb	WB,IF	H,M,R
A19631	ACTN3 Rabbit mAb	WB,IF	H,M,R
A2279	ACVR1B Rabbit mAb	WB,IHC,IF	H,M,R
A9755	ACVRL1 Rabbit mAb	WB,IF	H,M,R
A6813	ADAM15 Rabbit mAb	WB,IHC,IF	H,M,R
A3370	ADAMTS13 Rabbit mAb	WB,IF	H,M,R
A11466	ADAR1 Rabbit mAb	WB,IHC	H,M,R
A2383	Adenylate kinase 4 Rabbit mAb	WB,IHC,IF	H,R
A5219	ADORA1 Rabbit mAb	WB,IHC	H,M,R
A12411	AGR2 Rabbit mAb	WB,IF	H,M,R
A4000	AhR Rabbit mAb	WB,IHC	H
A19536	AIF Rabbit mAb	WB,IHC	H,M,R
A19776	AIF1/IBA1 Rabbit mAb	WB,IHC,IF	H,M,R
A3884	AKR1C3 Rabbit mAb	WB,IF	H,M
A17909	AKT1 Rabbit mAb	WB,IHC,IF	H,M,R
A2696	AKT1/2 Rabbit mAb	WB,IHC	H,M,R
A11625	Albumin Rabbit mAb	WB	H,M,R,Cow
A0157	ALDH1A1 Rabbit mAb	WB,IF,IP	H,M,R

货号	名称	应用	物种
A9123	ALDH1A2 Rabbit mAb	WB,IHC	H,M,R
A7707	ALDH1L1 Rabbit mAb	WB,IHC,IF	H,M,R
A11445	ALDOA Rabbit mAb	WB,IHC	H,M,R
A0514	Alkaline Phosphatase (ALPL) Rabbit mAb	WB,IHC	H,M,R
A9221	ALKBH1 Rabbit mAb	WB,IHC,IF	H,M,R
A2877	ALOX5 Rabbit mAb	WB,IHC,IF	H,M,R
A4899	alpha + beta Synuclein Rabbit mAb	WB,IF	H,M,R
A2418	Alpha Fodrin Rabbit mAb	WB,IHC	H,M,R
A3596	Alpha Internexin Rabbit mAb	WB,IHC,IF	H,M,R
A19754	alpha Sarcoglycan (SGCA) Rabbit mAb	WB,IHC,IF	H,M,R
A9752	Alpha-2-Macroglobulin (A2M) Rabbit mAb	WB,IHC,IF	H,M,R
A4344	AMPK β 1 Rabbit mAb	WB,IHC,IF	H,M,R
A19611	Androgen Receptor Rabbit mAb	WB,IHC,IF	H,M,R
A3757	ANGPT1 Rabbit mAb	WB,IHC	H,M,R
A11306	ANGPT2 Rabbit mAb	WB,IHC	H,M,R
A5225	ANGPTL3 Rabbit mAb	WB,IHC,IF	H,M,R
A11235	Annexin A2 Rabbit mAb	WB,IHC	H,M,R
A9203	Annexin A4/Annexin IV Rabbit mAb	WB,IF	H,M,R
A3489	ANP32B Rabbit mAb	WB,IHC,IF	H,M,R
A19269	AP1G1 Rabbit mAb	WB,IHC	H,M,R
A4403	AP2A1 Rabbit mAb	WB,IF	H,M,R
A11070	AP2M1 Rabbit mAb	WB,IF	H,M,R
A11238	APEX1/APE1 Rabbit mAb	WB,IHC	H,M,R
A4412	APOL1 Rabbit mAb	WB,IF	H,M
A4163	Apolipoprotein A1 Rabbit mAb	WB,IF	H
A3321	Apolipoprotein A2 Rabbit mAb	WB,IF	H,M,R
A2209	Apolipoprotein C3 Rabbit mAb	WB,IHC,IF	H,M,R
A4195	Aquaporin-1 (AQP1) Rabbit mAb	WB,IHC	H,M,R
A8687	A-RAF Rabbit mAb	WB,IHC	H
A9195	ARF1 Rabbit mAb	WB,IF	H,M,R
A3382	ARF5 Rabbit mAb	WB,IF	H,M,R
A11485	ARF6 Rabbit mAb	WB,IF	H,M,R
A4923	Arginase 1 (ARG1) Rabbit mAb	WB,IHC,IF	H,M,R
A19233	Arginase 2 (ARG2) Rabbit mAb	WB,IHC	H

货号	名称	应用	物种
A19709	Argonaute-2 Rabbit mAb	WB,IHC	H,M,R
A12238	Aromatase (CYP19A1) Rabbit mAb	WB,IHC,IF	H
A4514	Arp3 Rabbit mAb	WB,IHC,IF	H,M,R
A1382	ARPC2 Rabbit mAb	WB,IF,IP	H,M,R
A4892	ASH2L Rabbit mAb	WB,IHC,IF	H,M,R
A9084	ASS1 Rabbit mAb	WB,IF	H,M,R
A3563	ATF5 Rabbit mAb	WB,IHC	H,M,R
A6848	ATG10 Rabbit mAb	WB,IF	H,M
A19610	ATG12 Rabbit mAb	WB,IHC,IP	H,M
A19594	ATG3 Rabbit mAb	WB,IHC	H,M,R
A5059	ATG4B Rabbit mAb	WB,IF	H,M,R
A19604	ATG7 Rabbit mAb	WB,IHC,IP	H,M,R
A19650	ATM Rabbit mAb	WB,IHC	H,M,R
A5113	ATP1B1 Rabbit mAb	WB,IHC	H,M,R
A19637	ATP1B3 Rabbit mAb	WB,IHC	H,M,R
A19280	ATP4B Rabbit mAb	WB,IHC	H,M,R
A11217	ATP5A1 Rabbit mAb	WB,IHC,IF	H,M,R
A11214	ATPB Rabbit mAb	WB,IHC,IF	H,M,R
A19539	Aurora B Rabbit mAb	WB,IF	H,R
A19595	Bad Rabbit mAb	WB,IHC	H,M,R
A3814	BAF57/SMARCE1 Rabbit mAb	WB,IF	H,M,R
A8775	BAG2 Rabbit mAb	WB,IHC	H,M,R
A4720	BAG3 Rabbit mAb	WB,IHC	H,M
A3726	BANF1 Rabbit mAb	WB,IHC,IF	H,M,R
A2259	BAP31 Rabbit mAb	WB,IHC,IF	H,M,R
A19693	Bcl-2 Rabbit mAb	WB,IHC,IP	H,M
A19808	BCL2L12 Rabbit mAb	WB,IHC	H
A19292	BCL2L14/BCLG Rabbit mAb	WB,IHC	H,M,R
A11642	beta 2 Microglobulin Rabbit mAb	WB,IHC	H,M,R
A19702	Bim Rabbit mAb	WB	H,M,R
A20684	Biotin Rabbit mAb	WB	Other (Wide Range)
A4908	BiP/GRP78 Rabbit mAb	WB	H,M
A11405	BMP4 Rabbit mAb	WB,IHC,IF	H
A4538	BMP6 Rabbit mAb	WB,IF	H,M,R

货号	名称	应用	物种
A4405	BNIP1 Rabbit mAb	WB,IHC	H,M,R
A19593	BNIP3 Rabbit mAb	WB,IHC	H,M,R
A12677	BRD4 Rabbit mAb	WB,IHC	H
A19556	BRG1/SMARCA4 Rabbit mAb	WB,IHC,IF	H,M,R
A19002	BTK Rabbit mAb	WB,IHC,IF	H,M
A8831	BUB3 Rabbit mAb	WB,IHC,IF	H,M,R
A6849	CAD Rabbit mAb	WB,IHC,IF	H,M
A19644	Cadherin 16 Rabbit mAb	WB,IF	H,M
A3244	Cadherin-13 (CDH13) Rabbit mAb	WB,IHC	H,M,R
A4284	Calbindin Rabbit mAb	WB,IF	M,R
A4884	Caldesmon Rabbit mAb	WB,IHC	H,R
A19643	Calmeglin Rabbit mAb	WB,IF	H,M,R
A19648	Calmeglin Rabbit mAb	WB,IHC,IF	H,M,R
A4885	Calmodulin 1/2/3 Rabbit mAb	WB,IF	H,M,R
A8710	Calpain 1 Rabbit mAb	WB,IF	H,M,R
A3734	Calponin Rabbit mAb	WB,IHC,IF	H,M,R
A11563	Calreticulin Rabbit mAb	WB,IHC	H,M,R
A19640	Calsequestrin 1 Rabbit mAb	WB,IF	H,M,R
A9196	CaMKII delta Rabbit mAb	WB,IF	H,M,R
A0198	CaMKII Rabbit mAb	WB,IHC	H,M,R
A9271	CAMKIV Rabbit mAb	WB,IHC,IF	H,M,R
A9148	Carbonic Anhydrase 1 (CA2) Rabbit mAb	WB,IHC	H,M,R
A4934	Carbonic Anhydrase 9 (CA9/G250) Rabbit mAb	WB,IHC	H
A9652	CARD11 Rabbit mAb	WB,IHC,IF	H,M,R
A0438	CARS Rabbit mAb	WB,IHC	H,R
A9308	Casein Kinase 1 alpha (CSNK1A1) Rabbit mAb	WB,IHC	H,M,R
A4519	Casein Kinase 2 beta (CSNK2B) Rabbit mAb	WB,IHC,IF	H,M,R
A19664	Caspase-3 p12 Rabbit mAb	WB,IHC,IF	H,M,R
A19552	Caspase-6 Rabbit mAb	WB,IHC,IF	H,M,R
A19559	Caspase-6 Rabbit mAb	WB,IF	H,M
A4919	Caspase-9 Rabbit mAb	WB	H
A11220	Catalase Rabbit mAb	WB,IHC	H,M,R
A19005	Cathepsin B Rabbit mAb	WB,IHC	H
A19680	Cathepsin D Rabbit mAb	WB,IHC	H,M

货号	名称	应用	物种
A4986	Cathepsin L/V/K/H Rabbit mAb	WB,IHC	H,M,R
A3762	Cation-independent M6PR (IGF2R) Rabbit mAb	WB,IF	H,M,R
A19006	Caveolin-1 Rabbit mAb	WB,IHC	H,M,R
A4890	Caveolin-2 Rabbit mAb	WB,IHC	H,M
A4891	Caveolin-3 Rabbit mAb	WB,IHC	H,M,R
A19263	CBR1 Rabbit mAb	WB,IHC	H,M,R
A11612	CBS Rabbit mAb	WB,IHC,IF	H,M,R
A0121	CCR7 Rabbit mAb	WB,IHC,IF	H,M,R
A4288	CCR8 Rabbit mAb	WB,IF	H,M
A19008	CD105/Endoglin Rabbit mAb	WB,IHC	H
A19009	CD11a/LFA-1A/ITGAL Rabbit mAb	WB,IHC	H
A19010	CD11b/ITGAM Rabbit mAb	WB,IHC	H
A11678	CD127/IL7R Rabbit mAb	WB,IF	H,M,R
A11669	CD13/ANPEP Rabbit mAb	WB,IHC,IF	H,M,R
A0818	CD133 Rabbit mAb	WB,IHC	H
A4174	CD138/Syndecan-1 Rabbit mAb	WB,IHC,IF	H,M,R
A19011	CD14 Rabbit mAb	WB,IHC,IF	H,M
A9703	CD146/MCAM Rabbit mAb	WB,IHC	H,M,R
A9545	CD172a/SIRP α Rabbit mAb	WB,IF	H,M,R
A19013	CD19 Rabbit mAb	WB,IHC	H
A4893	CD20 Rabbit mAb	WB,IHC,IF	H
A9649	CD209 Rabbit mAb	WB,IHC,IF	H,M,R
A19724	CD239/BCAM Rabbit mAb	WB,IHC	H,M,R
A11505	CD27 Rabbit mAb	WB,IF	H,R
A9789	CD30 Rabbit mAb	WB,IF	H,R
A18643	CD31/PECAM1 Rabbit mAb	WB,IF,IP	H
A19014	CD31/PECAM1 Rabbit mAb	WB,IHC	H
A20228	CD31/PECAM1 Rabbit mAb	WB,IHC,IP	H
A4900	CD31/PECAM1 Rabbit mAb	WB,IF	H,M
A19015	CD34 Rabbit mAb	WB,IHC	H,M,R
A9770	CD3D Rabbit mAb	WB,IF	H,M,R
A19017	CD3E Rabbit mAb	WB,IHC	H
A4085	CD3G Rabbit mAb	WB,IHC,IF	H,R
A11157	CD3H Rabbit mAb	WB,IHC	H

货号	名称	应用	物种
A19018	CD4 Rabbit mAb	WB,IHC	H
A19019	CD40L Rabbit mAb	WB	M,R
A11490	CD41/ITGA2B Rabbit mAb	WB,IHC	H,M,R
A19021	CD45 Rabbit mAb	WB,IHC	H
A19022	CD46 Rabbit mAb	WB,IHC	H
A9557	CD5 Rabbit mAb	WB,IHC,IF	H
A11283	CD55 Rabbit mAb	WB,IHC	H
A4090	CD59 Rabbit mAb	WB,IHC	H,M,R
A4989	CD62P/P-selectin Rabbit mAb	WB,IHC	H,M,R
A9560	CD7 Rabbit mAb	WB,IHC	H
A9149	CD74 Rabbit mAb	WB,IHC	H
A19024	CD79a Rabbit mAb	WB,IHC	H
A9264	CD82 Rabbit mAb	WB,IHC,IF	H,M
A4234	CD83 Rabbit mAb	WB,IF	H,M,R
A0663	CD8A Rabbit mAb	WB,IHC	H
A19027	CD9 Rabbit mAb	WB,IHC	H,M
A12623	CD90/Thy1 Rabbit mAb	WB,IHC	H,R
A6025	CDC23 Rabbit mAb	WB,IF	H,M,R
A9758	CDC25B Rabbit mAb	WB,IHC	H
A12234	Cdc25C Rabbit mAb	WB,IF	H,M,R
A19794	CDC40 Rabbit mAb	WB,IHC	H,M,R
A19731	CDC42EP1 Rabbit mAb	WB,IHC	H,M,R
A9527	CDC5L Rabbit mAb	WB,IHC	H,M
A4921	CDCA3 Rabbit mAb	WB,IHC	H,M,R
A4044	CDCA5 Rabbit mAb	WB,IHC	H
A11420	CDK1 Rabbit mAb	WB,IHC,IF	H,M,R
A11145	CDK9 Rabbit mAb	WB,IHC,IF,IP,ChIP-seq	H,M,R
A8751	CDKN2C/p18-INK4C Rabbit mAb	WB,IF	H,M,R
A19030	CDX2 Rabbit mAb	WB,IHC,IF,IP	H,M,R
A11626	CEACAM1 Rabbit mAb	WB,IHC	H,R
A5971	CEACAM6 Rabbit mAb	WB,IF	H,M,R
A5930	CEND1 Rabbit mAb	WB,IF	H,M,R
A3975	CENPC Rabbit mAb	WB,IF	H,M,R
A20229	Ceruloplasmin Rabbit mAb	WB,IHC,IF	H,M,R

货号	名称	应用	物种
A11478	CES1 Rabbit mAb	WB,IF	H,M,R
A19031	CHAT Rabbit mAb	WB,IF	H
A6118	CHD3 Rabbit mAb	WB,IHC,IF	H,M,R
A4194	Chk1 Rabbit mAb	WB,IHC	H
A19543	Chk2 Rabbit mAb	WB,IHC	H
A19244	CHMP2B Rabbit mAb	WB,IHC	H,M,R
A9576	Chromogranin A Rabbit mAb	WB,IHC,IF	H,M,R
A19688	clAP1 Rabbit mAb	WB,IHC	H,R
A19290	CIRBP Rabbit mAb	WB,IHC,IF	H,M,R
A4569	CitRe synthetase Rabbit mAb	WB,IHC,IF	H,M,R
A4943	Clathrin heavy chain Rabbit mAb	WB,IHC,IF	H,M,R
A4387	CLEC3B/Tetranectin Rabbit mAb	WB,IF	H,M,R
A7265	CLOCK Rabbit mAb	WB,IF	H,M,R
A3214	CLPP Rabbit mAb	WB,IHC,IF	H,M,R
A19033	CNPase Rabbit mAb	WB,IHC,IF	H,M,R
A4532	Coagulation Protein C Rabbit mAb	WB,IF	H,M,R
A0817	Collagen III alpha 1/COL3A1 Rabbit mAb	WB,IF,IP	H,M,R
A9738	Collagen VI/COL6A1 Rabbit mAb	WB,IHC	H,M,R
A4808	Collagen XVII/COL17A1 Rabbit mAb	WB,IHC	H
A3545	Complement C4A Rabbit mAb	WB,IF	H
A19651	COPA Rabbit mAb	WB,IHC	H,M,R
A19584	COPS3 Rabbit mAb	WB,IHC	H,M,R
A3568	CoREST/RCOR1 Rabbit mAb	WB,IHC	H,M,R
A9518	Cortactin Rabbit mAb	WB,IHC,IF	H,M,R
A11631	COX IV Rabbit mAb	WB,IHC	H,M,R
A3560	COX2/PTGS2 Rabbit mAb	WB,IHC	H,M,R
A2228	COX6B1 Rabbit mAb	WB,IHC	H
A4214	CPS1 Rabbit mAb	WB,IHC,IF	H,M,R
A1267	CR2/CD21 Rabbit mAb	WB,IF	H,R
A5434	CRABP1 Rabbit mAb	WB,IHC,IF	H,M,R
A9265	CRALBP Rabbit mAb	WB,IHC,IF	H,M,R
A19003	C-Reactive Protein (CRP) Rabbit mAb	WB,IHC	H,M
A4907	Creatine kinase B type (CKB) Rabbit mAb	WB,IHC	H,M,R
A19625	CRM1/XPO1 Rabbit mAb	WB,IHC,IF	H,M,R

货号	名称	应用	物种
A4411	CRMP2/DPYSL2 Rabbit mAb	WB,IF	H,M,R
A3503	CRMP5/DPYSL5 Rabbit mAb	WB,IHC	H,M,R
A5111	CRYAA Rabbit mAb	WB,IF	M,R
A9633	CRYAB Rabbit mAb	WB,IHC,IF	H,R
A0666	CSK Rabbit mAb	WB,IHC	H,M,R
A0463	CTBP2 Rabbit mAb	IHC,IF	H,M,R
A9231	CTIP1/BCL11A Rabbit mAb	WB,IHC	H,M,R
A19804	CTIP2/BCL11B Rabbit mAb	WB,IHC	H,M
A19282	CTNNA3 Rabbit mAb	WB,IHC	H,M,R
A0773	CTR1/SLC31A1 Rabbit mAb	WB	H,M,R
A19660	CTRL Rabbit mAb	WB,IHC,IF	H,M,R
A5076	Cullin 2 Rabbit mAb	WB,IHC	H,M,R
A19623	Cullin 3 Rabbit mAb	WB,IHC	H,M,R
A5064	Cullin 4A Rabbit mAb	WB,IHC	H
A3694	CXCL4/PF4 Rabbit mAb	WB,IHC,IF	H,M,R
A11294	CXCR3 Rabbit mAb	WB,IF	H,M,R
A8950	CXCR5 Rabbit mAb	WB,IF	H,M,R
A19036	Cyclin A2 Rabbit mAb	WB,IHC,IF	H,M,R
A7956	Cyclin B2 Rabbit mAb	WB,IHC,IF	H,M,R
A4076	Cyclin H Rabbit mAb	WB,IHC,IF	H,M
A5097	Cyclophilin 40 Rabbit mAb	WB,IHC	H,M,R
A5067	CYP17A1 Rabbit mAb	WB,IF	H,M
A5982	CYP26A1 Rabbit mAb	WB,IHC,IF	H,M,R
A9562	CYP2D6 Rabbit mAb	WB,IF	H,M,R
A5032	CYPOR Rabbit mAb	WB,IHC,IF	H,M,R
A19662	Cytochrome P450 4A (CYP4A11) Rabbit mAb	WB,IHC	H,M,R
A9776	CytokeRin 1 (KRT1) Rabbit mAb	WB,IHC,IF	H,M,R
A4669	CytokeRin 10 Rabbit mAb	WB,IF	H,M,R
A0411	CytokeRin 13 (KRT13) Rabbit mAb	WB,IHC,IF	H,M,R
A19039	CytokeRin 14 (KRT14) Rabbit mAb	WB,IHC	H,M,R
A4854	CytokeRin 15 (KRT15) Rabbit mAb	WB,IF	H,R
A8690	CytokeRin 16 (KRT16) Rabbit mAb	WB,IHC,IF	H,M,R
A19778	CytokeRin 18 (KRT18) Rabbit mAb	WB,IHC,IF	H,M,R
A19041	CytokeRin 20 (KRT20) Rabbit mAb	WB,IHC	H,R

货号	名称	应用	物种
A2615	Cytokerin 2e (KRT2) Rabbit mAb	WB,IHC,IF	H,M,R
A0013	Cytokerin 4 (KRT4) Rabbit mAb	WB,IHC,IF	H,M,R
A11396	Cytokerin 5 (KRT5) Rabbit mAb	WB,IHC,IF	H,M,R
A19827	Cytokerin 6 (KRT6) Rabbit mAb	WB,IHC	H
A4631	Cytokerin 6 (KRT6) Rabbit mAb	WB,IHC	H,M,R
A4357	Cytokerin 7 (KRT7) Rabbit mAb	WB,IHC	H,M,R
A4765	Cytokerin 7 (KRT7) Rabbit mAb	WB,IHC	H,M,R
A6824	DCP1A Rabbit mAb	WB,IHC,IF	H,M
A3581	DCR2 Rabbit mAb	WB,IHC	H,M,R
A5022	DDB1 Rabbit mAb	WB,IHC	H,M,R
A5731	DDX17 Rabbit mAb	WB,IHC	H,M,R
A11339	DDX5 Rabbit mAb	WB,IHC,IF	H,M,R
A0315	DEK Rabbit mAb	WB,IHC,IF	H,R
A3736	Desmin Rabbit mAb	WB,IHC,IF	H,M,R
A5227	Desmoglein 3 Rabbit mAb	WB,IF	H,M
A4563	DHX9/RNA Helicase A Rabbit mAb	WB,IHC,IF,IP	H,M,R
A20267	Digoxin Rabbit mAb	DB	H,M,R,Other (Wide Range)
A4678	DISC1 Rabbit mAb	WB,IF	H,M,R
A4407	DKC1 Rabbit mAb	WB,IHC	H,M,R
A3892	DKK3 Rabbit mAb	WB,IF	H,M,R
A5220	DLDH/DLD Rabbit mAb	WB,IHC,IF	H,M,R
A2412	DNA Polymerase beta Rabbit mAb	WB,IHC	H,M,R
A12409	DNA topoisomerase I (TOP1) Rabbit mAb	WB,IHC,IF	H,M,R
A4825	DNAJC15 Rabbit mAb	WB,IF	H,M,R
A2382	DNPH1 Rabbit mAb	WB,IF	H,M
A19250	DOK1 Rabbit mAb	WB,IF	H,M,R
A4587	Dopamine Receptor D3 (DRD3) Rabbit mAb	WB	H,M,R
A5214	DP1/TFDP1 Rabbit mAb	WB,IHC	H,R
A19043	DR5/TRAILR2/TNFRSF10B Rabbit mAb	WB,IHC,IF	H,M,R
A19598	DROSHA Rabbit mAb	WB,IF	H
A0133	DUSP6 Rabbit mAb	WB,IHC	H,M,R
A19260	Dynamin 3 Rabbit mAb	WB	H,M
A4353	DYNLL1 Rabbit mAb	WB,IHC	H,M,R

货号	名称	应用	物种
A19572	Dysferlin (Romeo) Rabbit mAb	WB,IHC	H,M
A9712	EAAT1/SLC1A3 Rabbit mAb	WB,IHC	H,M,R
A19768	EB3/MAPRE3 Rabbit mAb	WB,IF	H,M,R
A19613	EBI3 Rabbit mAb	WB,IF	H,M,R
A5114	EBP1/PA2G4 Rabbit mAb	WB,IHC,IF	H,M,R
A11545	eEF1A1 Rabbit mAb	WB,IF	H,M,R
A8744	eEF1G Rabbit mAb	WB,IHC,IF	H,M,R
A9721	eEF2 Rabbit mAb	WB,IHC,IF	H,M,R
A4659	EG5/KIF11 Rabbit mAb	WB	H
A3496	EHD1 Rabbit mAb	WB,IHC	H,M,R
A9709	eIF2A Rabbit mAb	WB,IHC	H,M,R
A8670	eIF2B epsilon (EIF2B5) Rabbit mAb	WB,IHC,IF	H,M,R
A19819	eIF2C3 Rabbit mAb	WB,IHC	H,M,R
A9143	eIF3B Rabbit mAb	WB,IF	H,M,R
A3431	eIF3e Rabbit mAb	WB,IHC	H,M,R
A11584	eIF4A1 Rabbit mAb	WB,IF,IP	H,M,R
A8985	eIF4A3 Rabbit mAb	WB,IF	H,M,R
A19045	eIF4EBP1 Rabbit mAb	WB,IHC	H,R
A5169	ELMO1 Rabbit mAb	WB,IHC	H,M,R
A4187	Emerin/EMD Rabbit mAb	WB,IHC,IF	H,M,R
A19647	ENTPD5 Rabbit mAb	WB,IHC	H,M,R
A4149	Epac1/RAPGEF3 Rabbit mAb	WB,IHC	H,M,R
A1107	EpCAM Rabbit mAb	WB,IHC	H,M,R
A19301	EpCAM Rabbit mAb	WB,IHC,IF	H,M,R
A9132	Ephrin A1 Rabbit mAb	WB,IHC,IF	H,M
A11349	Ephrin B2 Rabbit mAb	WB,IHC	H,M,R
A5978	ERAP1 Rabbit mAb	WB,IHC	H,M,R
A19561	ERK1 Rabbit mAb	WB,IHC,IP	H,M,R
A4782	ERK1/2 Rabbit mAb	WB	H,M,R
A3948	ERK5 Rabbit mAb	WB,IHC	H,M,R
A0781	Erlin-2 Rabbit mAb	WB,IHC	H,M,R
A19771	ERp29 Rabbit mAb	WB,IHC	H,M,R
A4376	ERp57 Rabbit mAb	WB,IHC	H
A9640	EWSR1 Rabbit mAb	WB,IHC,IF	H,M,R

货号	名称	应用	物种
A0341	EXOSC7 Rabbit mAb	WB,IHC	H
A3813	Exportin 5 (XPO5) Rabbit mAb	WB,IHC	H,M,R
A19048	Ezrin Rabbit mAb	WB,IHC	H
A4099	FAAH1 Rabbit mAb	WB,IF	H,M,R
A11213	FABP1 Rabbit mAb	WB,IF	H,M,R
A8789	FABP3 Rabbit mAb	WB,IF	M,R
A8895	FABPI Rabbit mAb	WB,IHC,IF	M,R
A19838	FADD Rabbit mAb	WB,IHC	H,M,R
A9566	Fascin/FSCN1 Rabbit mAb	WB,IHC	H,M,R
A9005	FAT10/UBD Rabbit mAb	WB,IF	H,M,R
A19050	Fatty Acid Synthase (FASN) Rabbit mAb	WB,IHC,IF	H,M
A11664	FBP1 Rabbit mAb	WB,IHC,IF	H,M,R
A3699	Fbx32/FBOX32 Rabbit mAb	WB,IHC	M,R
A8999	FEN-1 Rabbit mAb	WB,IHC	H,M
A11241	Ferritin Light Chain Rabbit mAb	WB,IHC	H,M,R
A5900	FGF1 Rabbit mAb	WB,IF	H,M,R
A3908	FGF21 Rabbit mAb	WB,IHC	H,M,R
A19052	FGFR3 Rabbit mAb	WB,IHC	H,M,R
A9197	FGFR4 Rabbit mAb	WB,IHC	H,M
A3670	FHL2 Rabbit mAb	WB	H,M,R
A0850	Fibrillarin/U3 RNP Rabbit mAb	WB,IF	H,M,R
A19658	Fibrinogen alpha chain (FGA) Rabbit mAb	WB,IHC	H,M,R
A3738	Filamin A Rabbit mAb	WB,IHC	H,M,R
A8701	FKBP38/FKBP8 Rabbit mAb	WB,IF	H,M,R
A9090	FKBP51/FKBP5 Rabbit mAb	WB,IHC	H,M,R
A4356	FKBP52/FKBP4 Rabbit mAb	WB,IF	H,M,R
A3023	Flotillin 1 Rabbit mAb	WB,IHC	H,M,R
A19222	FMO3 Rabbit mAb	WB,IF	H,M,R
A4539	FMRP Rabbit mAb	WB,IHC	H,M,R
A2611	FNTB Rabbit mAb	WB,IHC	H,R
A19109	FSP1/S100A4 Rabbit mAb	WB,IHC,IF	H,M,R
A3861	FTO Rabbit mAb	WB,IHC,IF	H
A1385	Fukutin Rabbit mAb	WB,IHC	H,M,R
A8697	FXR1 Rabbit mAb	WB,IF	H,M,R

货号	名称	应用	物种
A3968	G3BP1 Rabbit mAb	WB,IF,IP	H,M,R
A9146	GABA transaminase (ABAT) Rabbit mAb	WB,IHC	H,M,R
A4732	Galectin 1/LGALS1 Rabbit mAb	WB,IHC,IF	H,M,R
A11198	Galectin 3/LGALS3 Rabbit mAb	WB,IHC	H,M,R
A3691	Galectin 4/LGALS4 Rabbit mAb	WB,IHC,IF	H,M,R
A2735	gamma-Synuclein Rabbit mAb	WB,IHC	H,M,R
A19055	GAP43 Rabbit mAb	WB,IF	H,M,R
A19056	GAPDH Rabbit mAb	WB,IHC,IF	H,M,R
AC041	GAPDH Rabbit mAb	WB,IHC	H,M,R
A0570	GBP1 Rabbit mAb	WB,IHC	H,R
A11292	GC1q R/C1QBP Rabbit mAb	WB,IHC	H,M,R
A4499	GCLC Rabbit mAb	WB,IHC	H,R
A2612	Gelsolin Rabbit mAb	WB,IHC	H,M
A19254	Gemin 2 Rabbit mAb	WB,IHC	H,M,R
A19058	GFAP Rabbit mAb	WB,IHC,IF	H,M,R
A19535	GFP Rabbit mAb	WB,IF	ALL
A3300	Gli3 Rabbit mAb	WB,IF	H,R
A19057	Glucosylceramidase beta (GBA) Rabbit mAb	WB,IHC	H,R
A5176	GLUD1 Rabbit mAb	WB,IHC,IF	H,M,R
A11316	GluR2/GRIA2 Rabbit mAb	WB,IHC,IF	H,M,R
A12307	GLUT2/SLC2A2 Rabbit mAb	WB,IF	H,M,R
A11043	Glutaminase (GLS) Rabbit mAb	WB,IF	H
A19641	Glutamine Synthetase (GLUL) Rabbit mAb	WB,IHC	H,M,R
A11557	Glutathione Synthetase (GSS) Rabbit mAb	WB,IF	H,M,R
A19221	Glycerol kinase (GK) Rabbit mAb	WB,IHC	H,M,R
A11472	Glycophorin C (GYPC) Rabbit mAb	WB,IHC	H,M,R
A11686	Glypican 3 (GPC3) Rabbit mAb	WB,IF	H
A9643	GNB2 Rabbit mAb	WB,IHC	H,M,R
A11363	GOT1 Rabbit mAb	WB,IHC,IF	H,M,R
A19245	GOT2 Rabbit mAb	WB,IHC	H,M,R
A19592	Granzyme B Rabbit mAb	WB	H
A2316	GRASP65 Rabbit mAb	WB,IF	H,M,R
A19059	GRB2 Rabbit mAb	WB,IHC,IF	H,M,R
A9722	GRIK2 Rabbit mAb	WB,IHC	H,M,R

货号	名称	应用	物种
A11731	GSK3 β Rabbit mAb	WB,IHC,IF	H,M,R
A19061	GST3/GSTP1 Rabbit mAb	WB,IHC	H,M,R
A9667	GSTK1 Rabbit mAb	WB,IF	H,M
A0536	GTF2F2 Rabbit mAb	WB,IHC	H,M,R
A4743	GTPBP4 Rabbit mAb	WB,IHC	H
A3687	Guanylyl Cyclase beta 1 (GUCY1B3) Rabbit mAb	WB,IHC,IF	H,M,R
A8778	Haptoglobin (HP) Rabbit mAb	WB,IHC,IF	H,M,R
A0670	HBG1 Rabbit mAb	WB,IF	H
A0704	HBG1 Rabbit mAb	WB,IHC,IF	H,M
A19626	HDAC2 Rabbit mAb	WB,IHC,IF,IP	H,M,R
A3572	HDAC6 Rabbit mAb	WB,IHC	H
A8865	HDAC8 Rabbit mAb	WB,IHC,IF	H,M,R
A2226	HDAC9 Rabbit mAb	WB,IHC,IF	H,R
A0589	HDGF Rabbit mAb	WB,IF	H,M,R
A3276	HE4/WFDC2 Rabbit mAb	WB,IF	H,M
A19062	Heme Oxygenase 1 (HO-1/HMOX1) Rabbit mAb	WB,IHC	H,M,R
A9293	Hemoglobin subunit alpha (HBA1) Rabbit mAb	WB,IHC,IF	H,M,R
A9133	Hemopexin (HPX) Rabbit mAb	WB,IHC,IF	H,M,R
A0925	HES1 Rabbit mAb	WB,IHC	M,R
A9768	HES5 Rabbit mAb	WB,IHC	H,M,R
A19532	HIF1 β /ARNT Rabbit mAb	WB,IHC	H
A9552	HIPK2 Rabbit mAb	WB,IF	H
A4342	Histone H1.0 Rabbit mAb	WB,IHC	H,M,R
A0646	Histone H1.2 Rabbit mAb	WB,IF	H,M,R
A3692	Histone H2A Rabbit mAb	WB,IHC	H,M,R,Other (Wide Range)
A11412	Histone H2AX Rabbit mAb	WB,IHC	H,M,R
A19812	Histone H2B Rabbit mAb	WB,IHC	H,M,R
A17562	Histone H3 Rabbit mAb	WB,IHC,IF	H,M,R
A4835	Histone H3.3 Rabbit mAb	WB,IHC	H,M,R,Other (Wide Range)
A19815	Histone H4 Rabbit mAb	WB,IHC	H,M,R
A11406	HLA-A Rabbit mAb	WB,IHC	H
A11252	HLA-DQA1 Rabbit mAb	WB,IHC	H,M,R
A10863	HLA-DRA Rabbit mAb	WB,IHC	H,R

货号	名称	应用	物种
A19835	HLA-DRB1 Rabbit mAb	WB,IHC	H
A4343	HMGA1 Rabbit mAb	WB,IHC,IF	H,M,R
A9168	HMGB2 Rabbit mAb	WB,IHC,IF	H,M,R
A19232	HMGCS2 Rabbit mAb	WB,IHC	H,M,R
A19053	HNF3 β /FOXA2 Rabbit mAb	WB,IF	H,R
A11564	hnRNP A1 Rabbit mAb	WB,IHC,IF,IP	H,M,R
A8958	hnRNP C Rabbit mAb	WB,IF	H,M
A19276	hnRNP E1/PCBP1 Rabbit mAb	WB	H,M,R
A0772	hnRNP K Rabbit mAb	WB,IHC,IF	H,M,R
A9609	hnRNP Q/SYNERIP Rabbit mAb	WB,IHC,IF	H,M,R
A4257	hnRNP U Rabbit mAb	WB,IHC,IF,IP	H,M,R
A19257	HOXA9 Rabbit mAb	WB,IHC	H,R
A3741	HP1 alpha/CBX5 Rabbit mAb	WB,IHC,IF	H,M,R
A4283	HP1 gamma/CBX3 Rabbit mAb	WB,IHC	H,M,R
A19619	HRAS Rabbit mAb	WB,IHC	H,M,R
A9091	HSP20/HSPB6 Rabbit mAb	WB,IHC,IF	H,M,R
A11156	HSP27/HSPB1 Rabbit mAb	WB,IF	H,M,R
A11698	HSP47/SERPINH1 Rabbit mAb	WB,IHC,IF	H,M
A0564	HSP60/HSPD1 Rabbit mAb	WB,IHC,IF	H,M,R
A4777	Hsp70 Rabbit mAb	WB,IF	H,M,R
A19574	Hsp90 β Rabbit mAb	WB,IF	H,M,R
A19825	HSPA12A Rabbit mAb	WB	H,M,R
A2286	HSPA14 Rabbit mAb	WB,IHC,IF	H,M
A6652	HSPA2 Rabbit mAb	WB,IF	H,M,R
A19705	HSPA4 Rabbit mAb	WB,IHC	H,M,R
A13518	HSPB8/HSP22 Rabbit mAb	WB,IF	H,M,R
A5580	HSPE1/HSP10/CPN10 Rabbit mAb	WB,IHC,IF	H,M,R
A19704	H IgA Rabbit mAb	WB	H
A19711	H IgG Rabbit mAb	WB	H,M
A19713	H IgG3 Rabbit mAb	WB	H,M,R
A19719	H IgM Rabbit mAb	WB,IHC	H
A19064	Huntingtin Rabbit mAb	WB,IHC,IF	H,M,R
A4827	Hydroxyl-Histone H2A-Y39 Rabbit mAb	WB,IHC	H,M,R
A19300	ICAM-1/CD54 Rabbit mAb	WB,IHC	H

货号	名称	应用	物种
A11270	IDH2 Rabbit mAb	WB,IHC,IF	H,M,R
A11653	IFNGR1 Rabbit mAb	WB,IF	H,M,R
A0830	IGF1 Rabbit mAb	WB,IHC	H,M,R
A5189	IGF2BP2/IMP2 Rabbit mAb	WB,IHC	H
A9018	IGFBP2 Rabbit mAb	WB,IHC	H
A4615	IGFBP7 Rabbit mAb	WB,IF	H,M,R
A3565	Ikaros Rabbit mAb	WB,IF	H,M
A19694	IKK α Rabbit mAb	WB	H,M
A19606	IKK β Rabbit mAb	WB,IHC	H
A3463	IKK ϵ Rabbit mAb	WB,IHC	H,M,R
A8614	IKZF3 Rabbit mAb	WB,IHC,IF	H,M,R
A12255	IL10 Rabbit mAb	WB,IF	H,M,R
A3524	IL13RA1 Rabbit mAb	WB,IHC	H
A11103	IL1RA Rabbit mAb	WB,IF	H,M,R
A19635	IL1 β Rabbit mAb	WB,IF	H
A9530	ILF2 Rabbit mAb	WB,IHC	H,M,R
A8186	ILF3 Rabbit mAb	WB,IHC	H,M,R
A3508	IMP3 Rabbit mAb	WB,IHC	H
A9208	IMPDH2 Rabbit mAb	WB,IF	H,M,R
A3774	iNOS Rabbit mAb	WB,IF	H,M
A19066	Insulin Rabbit mAb	WB,IHC,IF	M,R
A19068	Integrin alpha 2 (ITGA2/CD49b) Rabbit mAb	WB,IHC	H,M,R
A4054	Integrin alpha 4 (ITGA4/CD49d) Rabbit mAb	WB,IHC	H,M,R
A19069	Integrin alpha 5 (ITGA5/CD49e) Rabbit mAb	WB,IHC	H,M,R
A19071	Integrin alpha V (ITGAV/CD51) Rabbit mAb	WB,IHC	H,M,R
A19073	Integrin beta 3 (ITGB3/CD61) Rabbit mAb	WB,IHC	H,M,R
A4596	Integrin beta 4 (ITGB4/CD104) Rabbit mAb	WB	H,M
A19072	Integrin- β 1/CD29 Rabbit mAb	WB,IHC	H,M,R
A13311	Involucrin Rabbit mAb	WB,IHC	H,M,R
A4436	IP3R Rabbit mAb	WB,IHC,IF	M,R
A9765	IQGAP1 Rabbit mAb	WB,IHC	H,M
A4655	IRAK2 Rabbit mAb	WB,IHC	H,R
A11106	IRF5 Rabbit mAb	WB	H,M
A3209	IRF6 Rabbit mAb	WB,IF	H,M,R

货号	名称	应用	物种
A2416	ISG15 Rabbit mAb	WB,IF	H,M,R
A0871	Islet1 Rabbit mAb	WB,IHC	H,M,R
A19714	IκBα Rabbit mAb	WB,IHC	H,M,R
A4087	JAB1/CSN5/COPS5 Rabbit mAb	WB,IHC,IF	H,M,R
A11963	JAK1 Rabbit mAb	WB,IHC	H,M,R
A19629	JAK2 Rabbit mAb	WB	H,M
A5051	JNK1 Rabbit mAb	WB,IHC,IF	H,M,R
A4867	JNK1/2/3 Rabbit mAb	WB,IF	H,M,R
A4848	JunB Rabbit mAb	WB,IHC	H
A5496	JunD Rabbit mAb	WB,IHC	H,M,R
A3606	K48-linkage Specific Ubiquitin Rabbit mAb	WB,IF	H,M,R
A3770	Kaiso/ZBTB33 Rabbit mAb	WB,IF	H,M,R
A19568	KAP1/TRIM28 Rabbit mAb	WB,IHC,IP	H,M,R
A18740	Kappa light chain Rabbit mAb	WB,IHC,IF,IP	H
A4423	KAT1/HAT1 Rabbit mAb	WB,IHC,IF	H,M,R
A3390	KAT8/MYST1/MOF Rabbit mAb	WB,IHC,IF	H,M
A0020	KAT9/Elp3 Rabbit mAb	WB,IHC	H,M,R
A4689	K-Cadherin (CDH6) Rabbit mAb	WB	H,M,R
A8711	KDM1 / LSD1 Rabbit mAb	WB,IHC	H
A6670	KDM4B/JMJD2B Rabbit mAb	WB,IF	H,M,R
A3886	KHDRBS1/Sam68 Rabbit mAb	WB,IHC,IF	H,M,R
A9075	KHSRP Rabbit mAb	WB,IHC	H,M,R
A7370	KIF3A Rabbit mAb	WB,IHC,IF	H,M,R
A9080	KIF4A Rabbit mAb	WB,IHC	H,M,R
A5216	KIFAP3 Rabbit mAb	WB,IF	H,M,R
A0077	KIFC1 Rabbit mAb	WB,IHC,IF	H,M
A13673	KLF4 Rabbit mAb	WB,IHC,IF	H,M,R
A19779	KRAS+HRAS+NRAS Rabbit mAb	WB,IF,IP	H,M,R
A11223	Ku70 Rabbit mAb	WB,IHC,IF	H,M
A12338	Ku80 Rabbit mAb	WB,IHC,IF	H
A11495	Lamin B1 Rabbit mAb	WB,IHC	H,M,R
A5001	Lamin B2 Rabbit mAb	WB,IHC,IF	H,M,R
A4373	Laminin beta 1 Rabbit mAb	WB,IF	H,M,R
A0593	LAMP2 Rabbit mAb	WB,IHC	H,M,R

货号	名称	应用	物种
A0861	LDHA Rabbit mAb	WB,IF	H,M,R
A5131	LDHB Rabbit mAb	WB,IHC,IP	H,R
A2387	LEFTY1/2 Rabbit mAb	WB,IHC	H,M,R
A19076	Leptin Receptor Rabbit mAb	WB,IHC	H,M,R
A9128	LGR6 Rabbit mAb	WB,IF	H,R
A5286	LI Cadherin/Cadherin-17 Rabbit mAb	WB,IHC,IF	H,M,R
A9627	Lin28 Rabbit mAb	WB,IF	H,M,R
A9524	LIN28B Rabbit mAb	WB,IHC	H
A3696	LIS1 Rabbit mAb	WB,IF	H,M,R
A4941	LMAN1 Rabbit mAb	WB,IHC	H,M,R
A11504	LOX Rabbit mAb	WB,IHC	H,M,R
A4708	LOXL2 Rabbit mAb	WB,IHC	H
A4664	L-Plastin/LCP1 Rabbit mAb	WB,IHC	H
A19226	LPP Rabbit mAb	WB,IHC,IF	H,M,R
A19820	LSD2/KDM1B/AOF1 Rabbit mAb	WB	H,M,R
A3355	LSP1 Rabbit mAb	WB,IF	H,M,R
A8918	LTA4H Rabbit mAb	WB,IF	H,M,R
A11593	Lumican (LUM) Rabbit mAb	WB,IHC	H,M
A3974	LXR α Rabbit mAb	WB,IHC,IF	H,M,R
A0641	Lysozyme (LYZ) Rabbit mAb	WB,IF	H,M
A8648	LysRS/KARS Rabbit mAb	WB,IHC,IF	H,M,R
A9059	macroH2A.1 Rabbit mAb	WB,IHC,IF	H,M,R
A11469	MAD2/MAD2L1 Rabbit mAb	WB,IHC	H,M,R
A4630	MAD2B/MAD2L2 Rabbit mAb	WB,IHC,IF	H,M,R
A19727	Mammaglobin A Rabbit mAb	WB,IHC	H
A12319	MAP1LC3A Rabbit mAb	WB,IHC	H,M,R
A11480	Mast Cell Chymase (CMA1) Rabbit mAb	WB,IHC	H,M,R
A19801	Mast Cell Tryptase (TPSB2) Rabbit mAb	WB,IHC	H,M,R
A1027	Matrin 3 Rabbit mAb	WB,IHC,IF	H,M,R
A8905	MBD3 Rabbit mAb	WB,IHC	H,M,R
A5149	MBNL1 Rabbit mAb	WB,IF	H
A10993	MCL1 Rabbit mAb	WB,IHC,IF	H,M,R
A11475	MCM3 Rabbit mAb	WB,IHC,IF,IP	H,M,R
A5008	MCM5 Rabbit mAb	WB,IHC	H,M,R

货号	名称	应用	物种
A11325	MCM7 Rabbit mAb	WB,IHC,IF,IP	H,R
A0623	MED4 Rabbit mAb	WB,IHC	H,M,R
A2710	MEF2A+MEF2C Rabbit mAb	WB,IHC	H,M
A19565	MEK1 Rabbit mAb	WB,IHC,IF,IP	H,M,R
A3263	MEK3 Rabbit mAb	WB,IF	H,M,R
A12457	MEK5 Rabbit mAb	WB	H,M,R
A19739	MEKK2 Rabbit mAb	WB,IHC	H,M,R
A9180	Mesothelin Rabbit mAb	WB,IF	H,R
A3673	mGluR2 Rabbit mAb	WB,IF	M,R
A0880	MGST1 Rabbit mAb	WB,IHC,IF	H,M,R
A8754	MHC class I (HLA-A/B) Rabbit mAb	WB,IF	H
A11148	MMP13 Rabbit mAb	WB,IHC,IF	H
A11418	MMP3 Rabbit mAb	WB,IHC	H,M,R
A4915	Moesin Rabbit mAb	WB,IHC	H,M,R
A4105	Monoamine Oxidase A (MAOA) Rabbit mAb	WB,IHC	H,M,R
A11597	Monoamine Oxidase B (MAOB) Rabbit mAb	WB,IHC	H,M,R
A20680	MonoMethyl-Histone H3-K18 Rabbit pAb	WB,IHC,IF,IP,ChIP	H,M,R
A2358	MonoMethyl-Histone H3-K9 Rabbit pAb	WB,IHC,IF,IP,ChIP,ChIP-seq	H,M,R,Other (Wide Range)
A0452	MRPL28 Rabbit mAb	WB,IF,IP	H,M,R
A8740	MSH2 Rabbit mAb	WB,IHC	H,M,R
A8795	MSH6 Rabbit mAb	WB,IF	H,M,R
A2396	mSin3A Rabbit mAb	WB,IHC,IF	H,M,R
A4168	MSMB Rabbit mAb	WB,IF	H,M
A0109	MSP/MST1 Rabbit mAb	WB,IHC,IF	H,M,R
A2401	MSR1 Rabbit mAb	WB,IHC	H,M,R
A3843	MTCO2 Rabbit mAb	WB,IHC	H
A19081	MUC1 Rabbit mAb	WB,IHC	H,M,R
A4767	MUC2 Rabbit mAb	WB,IHC	H,M,R
A3438	MUC4 Rabbit mAb	WB,IHC	H,R
A9122	Musashi-1 (MSI1) Rabbit mAb	WB,IF	H,M
A19814	Musashi-2 (MSI2) Rabbit mAb	WB,IHC	H,M,R
A5154	Muscarinic AChR M2 (ACM2) Rabbit mAb	WB,IHC	H,M,R
AE070	Myc-Tag Rabbit mAb	WB,IF,FC,IP	

货号	名称	应用	物种
A11162	Myelin Basic Protein Rabbit mAb	WB,IF	M,R
A3992	Myelin oligodendrocyte glycoprotein Rabbit mAb	WB,IF	M,R
A8742	MYL2 Rabbit mAb	WB,IF	M,R
A4963	Myosin heavy chain Rabbit mAb	WB,IHC,IF	M,R
A19841	N6-methyladenosine / m6A Rabbit mAb	Nucleotide Array,DB,IF	H,M,R,Other (Wide Range)
A11683	Na ⁺ /K ⁺ -ATPase Rabbit mAb	WB,IHC,IF	H,M,R
A6174	NAP1L1 Rabbit mAb	WB,IF	H,M,R
A19286	NAT10 Rabbit mAb	WB,IHC	H,M,R
A0935	NCF4/p40-phox Rabbit mAb	WB,IHC,IP	H,M
A9129	NCK1 Rabbit mAb	WB,IHC,IF	H,M,R
A4050	NDRG1 Rabbit mAb	WB,IHC	H,M,R
A4007	NDUFAF1 Rabbit mAb	WB,IF	H
A19732	NDUFB8 Rabbit mAb	WB,IHC,IF	H,M,R
A4602	NDUFS3 Rabbit mAb	WB,IHC	H,M,R
A8691	NDUFS4 Rabbit mAb	WB,IHC,IF	H,M,R
A9622	Nectin 2/CD112 Rabbit mAb	WB,IHC	H,M,R
A3536	NEK6 Rabbit mAb	WB,IHC	H,M,R
A19781	Nesprin 1 Rabbit mAb	WB,IF	H,M,R
A19086	NeuN Rabbit mAb	WB,IHC,IF	H,M,R
A19084	Neurofilament H Rabbit mAb	WB,IHC,IF	H,M,R
A20269	Neurofilament L Rabbit mAb	WB,IHC,IF,IP	H,M,R
A19085	Neurofilament M Rabbit mAb	WB,IF	H
A8953	Neutrophil Elastase (ELANE) Rabbit mAb	WB,IHC,IF	H
A19605	NF-κB2 Rabbit mAb	WB,IHC	H,M,R
A19088	NGFR Rabbit mAb	WB,IHC	H,M,R
A19270	NLK Rabbit mAb	WB,IHC	H,M,R
A11699	NMDAR1 Rabbit mAb	WB,IF	H,M,R
A19089	NMDAR2A Rabbit mAb	WB,IHC	H,M,R
A2649	nNOS Rabbit mAb	WB,IF	H,M,R
A3800	NONO/p54nrb Rabbit mAb	WB,IHC,IF	H,M,R
A3703	NOXA2/p67phox Rabbit mAb	WB,IHC	M,R
A19586	NQO1 Rabbit mAb	WB,IHC,IP	H,M,R
A4735	NRAS Rabbit mAb	WB,IF	H,M,R

货号	名称	应用	物种
A3252	NRF1 Rabbit mAb	WB,IHC	H,M,R
A0609	NUDT5 Rabbit mAb	WB,IF	H,M,R
A4416	NUMA1 Rabbit mAb	WB,IHC	H,M,R
A2472	NUP153 Rabbit mAb	WB,IHC,IF	H,M,R
A2270	N-WASP/WASL Rabbit mAb	WB,IHC	H,M,R
A4997	OGG1 Rabbit mAb	WB,IHC	H,M,R
A8705	Oncostatin M Rabbit mAb	WB,IF	H,M,R
A19736	ORM1 Rabbit mAb	WB,IHC	H
A1042	ORP150/GRP170/HSP12A Rabbit mAb	WB,IHC	H,M,R
A19093	P Glycoprotein Rabbit mAb	WB,IHC	H,M,R
A5194	P23/PTGES3 Rabbit mAb	WB,IF	H,M,R
A5049	p38 MAPK Rabbit mAb	WB,IHC,IF	H,M,R
A19239	P4HB Rabbit mAb	WB,IHC	H,M,R
A5952	p53 DINP1 Rabbit mAb	WB,IHC,IF	H,M,R
A5073	p60 CAF1 Rabbit mAb	WB,IHC	H
A19652	p63 Rabbit mAb	WB,IHC	H,M
A9100	p70 S6 Kinase 2 Rabbit mAb	WB,IHC,IF	H,M,R
A2670	p73 Rabbit mAb	WB,IHC	H,M,R
A1735	PABPN1 Rabbit mAb	WB,IHC,IF	H,M,R
A3437	PAF1 Rabbit mAb	WB,IHC	H,M,R
A19608	PAK1 Rabbit mAb	WB,IF	H,M,R
A4553	PAK2 Rabbit mAb	WB,IF	H,M,R
A3363	PAK3 Rabbit mAb	WB,IHC	H,M,R
A4903	Pan Cadherin Rabbit mAb	WB,IHC	H,M,R
AP1162	Pan Phospho-Tyrosine Rabbit mAb	WB	H,M,R,Other (Wide Range)
A18675	Pan-AKT Rabbit mAb	WB,IHC,IF,IP	H,M,R
A5167	Pannexin 1 Rabbit mAb	WB,IHC	H,M,R
A5103	PAR2 Rabbit mAb	WB,IHC	H,M,R
A4983	PAR4 Rabbit mAb	WB,IHC,IF	H,M
A9704	PaRhyroid Hormone (PTH) Rabbit mAb	WB,IF	H,M,R
A19098	Parvalbumin (PVALB) Rabbit mAb	WB,IHC	H,M,R
A9607	PAX5 Rabbit mAb	WB,IHC	H
A19099	PAX6 Rabbit mAb	WB,IHC	H,M,R

货号	名称	应用	物种
A19741	PAX9 Rabbit mAb	WB,IHC	H,M,R
A4881	PBR/TSPO Rabbit mAb	WB,IHC,IF	H,M
A0334	PBRM1 Rabbit mAb	WB,IHC	H,M,R
A12427	PCNA Rabbit mAb	WB,IHC,IP	H,M,R
A9068	PDCD4 Rabbit mAb	WB,IHC	H,M,R
A1510	PDCD7 Rabbit mAb	WB,IHC	H,M,R
A19531	PDGFRB Rabbit mAb	WB,IF	H,M
A13687	PDHA1 Rabbit mAb	WB,IF	H,M,R
A4645	PDHB Rabbit mAb	WB,IHC	H,M,R
A4260	PDIA6 Rabbit mAb	WB,IF	H,M,R
A4737	PDK2 Rabbit mAb	WB,IHC	H,M,R
A19224	PDLIM1/CLP36 Rabbit mAb	WB,IF	H,M,R
A8678	PEN2/PSENNEN Rabbit mAb	WB,IF	H,M,R
A4466	PEPCK/PCK2 Rabbit mAb	WB,IF	H,M,R
A5107	PER2 Rabbit mAb	WB,IF	H,M
A2219	PER3 Rabbit mAb	WB,IF	H,M,R
A4758	Perilipin A Rabbit mAb	WB,IHC,IF	M,R
A2398	Peroxiredoxin 3 (PRDX3) Rabbit mAb	WB,IHC,IF	H,M
A9131	Peroxiredoxin 4 (PRDX4) Rabbit mAb	WB,IHC,IF	H,M,R
A5937	PERP Rabbit mAb	WB,IHC,IF	H,M,R
A4101	PFDN5 Rabbit mAb	WB,IF	H,M,R
A3934	PFKFB3 Rabbit mAb	WB,IHC	H,M,R
A4170	PGAM1 Rabbit mAb	WB,IHC	H,M,R
A0563	PGD Rabbit mAb	WB,IF,IP	H,M,R
A9053	PGE receptor EP2 (PTGER2) Rabbit mAb	WB,IHC,IF	H,M,R
A19101	PGP9.5/Uchl1 Rabbit mAb	WB,IHC,IF	H,M,R
A20380	PGP9.5/Uchl1 Rabbit mAb	WB,IHC	H,M,R
AP1259	Phospho-Akt (Thr308) Rabbit mAb	WB	H,M,R
AP0982	Phospho-AKT1-S124 Rabbit mAb	WB	H
AP0981	Phospho-AKT1-S129 Rabbit mAb	WB	H
AP0980	Phospho-AKT1-T450 Rabbit mAb	WB	H,R
AP1208	Phospho-Akt-S473 Rabbit mAb	WB	H,M,R
AP1002	Phospho-AMPK α 1-S496 Rabbit mAb	WB	H,M
AP1051	Phospho-ATF2-T71 Rabbit mAb	WB	H,M

货号	名称	应用	物种
AP1030	Phospho-ATM-S1981 Rabbit mAb	WB	H
AP1010	Phospho-Bad-S112 Rabbit mAb	WB	H,M
AP1041	Phospho-B-Raf-T401 Rabbit mAb	WB	H,M,R
AP1141	Phospho-CDC37-S13 Rabbit mAb	WB	H,M,R
AP1153	Phospho-Cdc6-S54 Rabbit mAb	WB	H,M,
AP1138	Phospho-CDK1/2/3-T14 Rabbit mAb	WB	H,M,R
AP1001	Phospho-Cdk1/2-T14 Rabbit mAb	WB	H,M,R
AP1005	Phospho-CDK2-Y15 Rabbit mAb	WB	H,M
AP1047	Phospho-Chk1-S296 Rabbit mAb	WB	H
AP0105	Phospho-c-Jun-S63 Rabbit mAb	WB,IF	H,M,R
AP1003	Phospho-c-Jun-T91 Rabbit mAb	WB,IHC	H,M,R
AP0989	Phospho-c-Myc-S62 Rabbit mAb	WB	H,M
AP0988	Phospho-c-Myc-T58/S62 Rabbit mAb	WB	H,R
AP1014	Phospho-Cyclin E1-T77 Rabbit mAb	WB	H,M
AP0994	Phospho-EGFR-Y1068 Rabbit mAb	WB	H
AP0992	Phospho-EGFR-Y1173 Rabbit mAb	WB	H
AP1052	Phospho-ErbB3/HER3-Y1222 Rabbit mAb	WB	H
AP0485	Phospho-ERK1-T202 + ERK2-T185 Rabbit mAb	WB,IHC	H,M,R
AP0974	Phospho-Erk1-T202/Y204 + Erk2-T185/Y187 Rabbit mAb	WB,IHC,IF	H,M,R
AP1004	Phospho-ESR α -S118 Rabbit mAb	WB	H,M
AP1088	Phospho-GSK3 β -S9 Rabbit mAb	WB,IHC	H,M
AP1129	Phospho-GYS1-S641 Rabbit mAb	WB	H,M
AP1132	Phospho-Histone H1.3-T17/Histone H1.4-T17 Rabbit mAb	WB	H,M,R
AP0687	Phospho-Histone H2AX-S139 Rabbit mAb	WB,IHC,IF	H,M,R
AP1152	Phospho-Histone H3.3-T3 Rabbit mAb	WB	H,M,R,Other (Wide Range)
AP0002	Phospho-Histone H3-S10 Rabbit mAb	WB,IHC,IF	H,M,R,Other (Wide Range)
AP1140	Phospho-HSF1-S326 Rabbit mAb	WB	H
AP1151	Phospho-HSL-S853 Rabbit mAb	WB	H,M
AP1032	Phospho-HSP27/HSPB1-S78 Rabbit mAb	WB	H,M
AP1031	Phospho-HSP27/HSPB1-S82 Rabbit mAb	WB	H,M,R
AP0995	Phospho-IRF3-S386 Rabbit mAb	WB	H
AP0707	Phospho-IkB α -S32 Rabbit mAb	WB	H,M,R

货号	名称	应用	物种
AP0999	Phospho-IkB α -S36 Rabbit mAb	WB,IHC	H,M,R
AP0631	Phospho-JNK1/2/3-T183/T183/T221 Rabbit mAb	WB,IHC,IF	H,M,R
AP1142	Phospho-JunD-S255 Rabbit mAb	WB	H,M
A3430	Phospholipid Scramblase 1 (PLSCR1) Rabbit mAb	WB,IHC,IF	H,M,R
AP1154	Phospho-LRRK2-S935 Rabbit mAb	WB	H,M,R
AP1050	Phospho-Lyn-Y396 Rabbit mAb	WB	H,M,R
AP1020	Phospho-MEK1-S298 Rabbit mAb	WB	H,R
AP1021	Phospho-MEK1-T292 Rabbit mAb	WB	H,M,R
AP0978	Phospho-mTOR-S2481 Rabbit mAb	WB	H,M
AP1133	Phospho-NRF2-S40 Rabbit mAb	WB,IHC	H,M,R
AP0984	Phospho-p53-S33 Rabbit mAb	WB	H
AP0987	Phospho-p53-S376 Rabbit mAb	WB	H,R
AP1137	Phospho-p53-S392 Rabbit mAb	WB	H,R
AP0985	Phospho-p53-S9 Rabbit mAb	WB	H,M
AP0986	Phospho-p53-T55 Rabbit mAb	WB	H
AP0502	Phospho-p70 S6 Kinase 1-T421/S424 Rabbit mAb	WB,IHC,IF	H,M,R
AP1147	Phospho-p90Rsk/RSK1/RPS6KA1-S380 Rabbit mAb	WB	R
AP1158	Phospho-PAK1/2/3-S144/S141/S139 Rabbit mAb	WB	H,M,R
AP1049	Phospho-PAK4/PAK5/PAK6-S474/S560/S602 Rabbit mAb	WB	H,M,R
AP1156	Phospho-Paxillin-Y118 Rabbit mAb	WB	H,M,R
AP1022	Phospho-PDHA1-S293 Rabbit mAb	WB,IF	H,M,R
AP1034	Phospho-PKA RI α (PRKAR2A)-S99 Rabbit mAb	WB	H,R
AP1045	Phospho-PKC alpha-T638 Rabbit mAb	WB	H,M,R
AP1149	Phospho-PKC zeta-T560 Rabbit mAb	WB	H,R
AP1134	Phospho-PKR/EIF2AK2-T446 Rabbit mAb	WB	H
AP1025	Phospho-PLK1-T210 Rabbit mAb	WB	H
AP0996	Phospho-POLR2A-S2 Rabbit mAb	WB	H,M,R
AP0997	Phospho-POLR2A-S5 Rabbit mAb	WB	H,M
AP1043	Phospho-PP2A Catalytic α -Y307 Rabbit mAb	WB	H,M
AP0998	Phospho-PTEN-S380 Rabbit mAb	WB	H
AP1018	Phospho-Rac1/Cdc42-S71 Rabbit mAb	WB	H
AP1011	Phospho-Raf1-S621 Rabbit mAb	WB	H,R
AP1040	Phospho-RPA2-T21 Rabbit mAb	WB	H,R

货号	名称	应用	物种
AP1007	Phospho-Smad2-S250 Rabbit mAb	WB	H,M
AP0727	Phospho-Smad3-S423/S425 Rabbit mAb	WB,IHC	H,M,R
AP1023	Phospho-Smad5-S463/S465 Rabbit mAb	WB	H
AP1145	Phospho-SMC1-S957 Rabbit mAb	WB,IHC	H,M,R
AP1027	Phospho-Src-Y419 Rabbit mAb	WB	H,M
AP1000	Phospho-STAT1-S727 Rabbit mAb	WB,IHC	H,M,R
AP0054	Phospho-STAT1-Y701 Rabbit mAb	WB	H,M,R
AP0715	Phospho-STAT3-S727 Rabbit mAb	WB,IHC	H,M,R
AP0705	Phospho-STAT3-Y705 Rabbit mAb	WB,IP	H,M
AP0758	Phospho-STAT5A-Y694 Rabbit mAb	WB,IHC	H,M,R
AP1028	Phospho-Tau-S396 Rabbit mAb	WB,IHC	H,M,R
AP0053	Phospho-Tau-T231 Rabbit mAb	WB,IF	H,M,R
AP1026	Phospho-TBK1/NAK-S172 Rabbit mAb	WB	H,M
AP1207	Phospho-VANGL2-S79/S82/S84 Rabbit mAb	WB	H,M,R
AP1206	Phospho-VANGL2-T78/S79/S82 Rabbit mAb	WB	H,M,R
A19742	PI3 Kinase p110 delta Rabbit mAb	WB,IHC	H,M,R
A4992	PI3 Kinase p85 alpha Rabbit mAb	WB,IHC	H,M,R
A4860	PI3 Kinase p85 beta Rabbit mAb	WB,IHC	H,M,R
A9670	PIAS1+PIAS2 Rabbit mAb	WB,IHC	H,M,R
A12295	PIK3C3/VPS34 Rabbit mAb	WB,IF	H,M,R
A5922	PIK3R4/VPS15 Rabbit mAb	WB,IF	H,M,R
A19695	PIM1 Rabbit mAb	WB,IHC	H,M,R
A0582	PIST/GOPC Rabbit mAb	WB,IHC,IF	H,M,R
A3490	PIWIL1 Rabbit mAb	WB,IHC,IF	M,R
A11107	PKC alpha Rabbit mAb	WB,IHC	H,M,R
A7778	PKC delta Rabbit mAb	WB,IHC,IF	H,M,R
A4304	Placental alkaline phosphatase (PLAP) Rabbit mAb	WB,IHC	H,M,R
A9644	PLTP Rabbit mAb	WB,IHC	H,M,R
A4172	PMP70/ABCD3 Rabbit mAb	WB,IF	H,M,R
A4577	PMS2 Rabbit mAb	WB,IHC	H
A4554	PODXL Rabbit mAb	WB,IHC	H
A3441	PON1 Rabbit mAb	WB,IHC,IF	H,M,R
A3683	POU2F1/OCT1 Rabbit mAb	WB,IHC	H,M
A4364	PP1 beta Rabbit mAb	WB,IHC,IF	H,M,R

货号	名称	应用	物种
A6175	PP2A alpha + beta Rabbit mAb	WB,IHC,IF	H,M,R
A19746	PRK2/PKN2 Rabbit mAb	WB,IF	H,M,R
A19745	PRKAR2B Rabbit mAb	WB,IHC	H,M,R
A8779	PRKRA/PACT Rabbit mAb	WB,IHC	H,M,R
A4502	PRMT1 Rabbit mAb	WB,IHC,IF	H,M,R
A19533	PRMT5 Rabbit mAb	WB,IHC	H,M,R
A19697	Progesterone Receptor Rabbit mAb	WB	H
A9144	Prohibitin 2 (PHB2) Rabbit mAb	WB,IHC,IF	H,M,R
A19530	Prohibitin Rabbit mAb	WB,IHC	H,M,R
A19268	PROK1 Rabbit mAb	WB,IF	H,M,R
A9047	PROX1 Rabbit mAb	WB,IHC	H,M,R
A9660	PRP19 Rabbit mAb	WB,IHC,IF	H,M,R
A4575	PRPF8 Rabbit mAb	WB,IHC,IF	H,M,R
A19275	PRSS2 Rabbit mAb	WB,IHC	H,M,R
A19669	PSD93/chapsyn-110/DLG2 Rabbit mAb	WB,IF	H,M,R
A9547	PSMA/FOLH1 Rabbit mAb	WB,IF	R
A3460	PSMA1 Rabbit mAb	WB,IF	H,M,R
A9182	PSMA2 Rabbit mAb	WB,IHC,IF	H,M,R
A5094	PSMA3 Rabbit mAb	WB,IHC	H,M,R
A9549	PSMB9/LMP2 Rabbit mAb	WB,IHC	H,M,R
A9608	PSMD14 Rabbit mAb	WB,IHC,IF	H,M,R
A3663	PSMD4 Rabbit mAb	WB,IF	H,M,R
A3487	PTBP1 Rabbit mAb	WB,IHC,IF	H,M,R
A9124	PTBP2 Rabbit mAb	WB,IHC,IF	H,M,R
A3752	PUMA Rabbit mAb	WB,IHC,IF	H,M,R
A0571	Pumilio 1 Rabbit mAb	WB,IF	H,M
A19783	Pumilio 2 Rabbit mAb	WB	H,M,R
A0193	QK1 Rabbit mAb	WB,IF,IP	H,M,R
A3251	RAB11A Rabbit mAb	WB,IF	H,M,R
A12308	RAB7 Rabbit mAb	WB,IF	H,M,R
A4894	RAB9A Rabbit mAb	WB,IHC,IF	H,M,R
AE078	Rabbit anti GFP-Tag mAb	WB,IF	
AE088	Rabbit anti GFP-Tag mAb	WB,IF	
AE077	Rabbit anti GST-Tag mAb	WB,IP	

货号	名称	应用	物种
AE086	Rabbit anti His-tag mAb	WB,IP	
AE089	Rabbit anti V5-Tag pAb	WB,IF,IP	
A5147	Rad23A Rabbit mAb	WB,IF	H,M,R
A5186	Rad52 Rabbit mAb	WB	H
A19638	Raf1 Rabbit mAb	WB	H,M,R
A8663	RAIDD/CRADD Rabbit mAb	WB,IF	H,R
A8664	RAIDD/CRADD Rabbit mAb	WB,IHC	H,R
A19238	RanBP9 Rabbit mAb	WB,IHC	H,M,R
A8758	RanGAP1 Rabbit mAb	WB,IHC,IF	H,M,R
A9725	RAP1A + RAP1B Rabbit mAb	WB,IF	H,M,R
A3621	Rap1GAP Rabbit mAb	WB,IF	H,M,R
A8992	Raptor Rabbit mAb	WB,IHC,IF	H,M,R
A19551	RAR α Rabbit mAb	WB,IHC	H,M,R
A4535	RAR β Rabbit mAb	WB,IF	H,M
A8807	RBP4 Rabbit mAb	WB,IF	H,M,R
A4081	RBPJK Rabbit mAb	WB,IHC	H,M
A0662	RCC1 Rabbit mAb	WB,IHC	H,M,R
A4855	RHOA/RHOB/RHOC Rabbit mAb	WB,IHC,IF	H,M,R
A7245	Rhodopsin Rabbit mAb	WB,IHC	M,R
A19750	RING1A Rabbit mAb	WB,IHC	H
A9598	RNF40 Rabbit mAb	WB,IHC	H,M
A2395	ROCK2 Rabbit mAb	WB,IHC	H,M,R
A9615	RPA70 Rabbit mAb	WB,IHC,IF	H,M,R
A9744	RPB3/POLR2C Rabbit mAb	WB,IHC	H,M,R
A3675	RPS19 Rabbit mAb	WB,IF	H,M,R
A9551	RPS20 Rabbit mAb	WB,IHC,IF	H,M,R
A4872	RPS3 Rabbit mAb	WB,IF	H,M,R
A5180	RUVBL1/TIP49A/PONTIN Rabbit mAb	WB,IHC	H,M,R
A19105	RXR α Rabbit mAb	WB,IF	H,M,R
A13614	S100A10 Rabbit mAb	WB,IHC	H
A3461	S100A6 Rabbit mAb	WB,IHC,IF	H,M,R
A12018	S100A8 Rabbit mAb	WB,IHC	H,M,R
A19108	S100B Rabbit mAb	WB,IHC,IF	H,M,R
A11874	S6 Ribosomal Protein (RPS6) Rabbit mAb	WB	H,M,R

货号	名称	应用	物种
A4607	SAMHD1 Rabbit mAb	WB,IHC	H
A19837	SATB2 Rabbit mAb	WB,IHC	H,M,R
A3491	SATB2 Rabbit mAb	WB,IHC,IF	H,M,R
A0796	SBDS Rabbit mAb	WB	H,M,R
A8613	Sec23A Rabbit mAb	WB,IHC,IF	H,M,R
A19615	Secretagogin (SCGN) Rabbit mAb	WB,IHC	H,M,R
A19639	SERCA1/ATP2A1 Rabbit mAb	WB,IF	H,M,R
A19110	Serotonin transporter Rabbit mAb	WB,IF	H,M,R
A11249	SERPINC1 Rabbit mAb	WB,IHC	H,M,R
A19787	sestrin 1 Rabbit mAb	WB,IHC	H,M,R
A9173	SET/TAF1 Rabbit mAb	WB,IHC,IF,IP	H,M,R
A9624	SF3B3/SAP130 Rabbit mAb	WB,IHC,IF	H,M,R
A4189	SFRP4 Rabbit mAb	WB,IHC,IF	H,M,R
A3571	SHIP1 Rabbit mAb	WB,IHC,IF	H,M,R
A19111	SHP1 Rabbit mAb	WB,IHC,IF	H,M,R
A19112	SHP2 Rabbit mAb	WB,IF,IP	H,M,R
A9021	SKP1 Rabbit mAb	WB,IF	H,M,R
A4046	SKP2 Rabbit mAb	WB,IHC	H,M,R
A19758	SLC22A3/OCT3 Rabbit mAb	WB,IHC	H,M
A3467	Slit2 Rabbit mAb	WB,IHC	H,M,R
A8889	Smac/Diablo Rabbit mAb	WB,IHC,IF	H,M,R
A3539	SMARCA5/SNF2H Rabbit mAb	WB,IF	H,M,R
A4275	SMARCC1/BAF155 Rabbit mAb	WB,IF	H,M
A4693	SMC1 Rabbit mAb	WB,IF	H,M,R
A19591	SMC3 Rabbit mAb	WB,IF	H,M,R
A3559	SMC4 Rabbit mAb	WB,IF	H,M,R
A4064	SMMHC/MYH11 Rabbit mAb	WB,IHC,IF	H,M,R
A9707	SMN1 Rabbit mAb	WB,IHC	H,R
A2278	SMURF2 Rabbit mAb	WB,IHC,IF,IP	H,M
A2262	SMYD3 Rabbit mAb	WB,IF	H,M,R
A0647	SNRPA1 Rabbit mAb	WB,IHC	H,M,R
A9190	SOCS2 Rabbit mAb	WB,IHC,IF	H,M,R
A19576	SOD2 Rabbit mAb	WB,IHC,IF	H,M,R
A12695	Sonic Hedgehog (Shh) Rabbit mAb	WB,IF	H,M,R

货号	名称	应用	物种
A8655	SOX10 Rabbit mAb	WB,IF	H,M,R
A8658	SOX10 Rabbit mAb	WB,IF	M,R
A19710	SOX9 Rabbit mAb	WB,IHC,IF,IP	H,M,R
A19649	SP1 Rabbit mAb	WB,IHC,IF	H,M,R
A9185	SR-B2/LIMPII Rabbit mAb	WB,IHC	H,M,R
A0827	SR-BI Rabbit mAb	WB,IHC	H,M,R
A19762	SRD5A2 Rabbit mAb	WB,IF	H,M,R
A4091	SRSF1/SF2/ASF Rabbit mAb	WB,IHC	H,M,R
A9054	SRSF3 Rabbit mAb	WB,IHC	H,M,R
A8698	SSTR5 Rabbit mAb	WB,IHC	H,M,R
A3588	STAT2 Rabbit mAb	WB	H,M,R
A4523	STAT4 Rabbit mAb	WB,IF	H,M,R
A19120	STAT6 Rabbit mAb	WB,IF	H,M
A4379	Stathmin 1 Rabbit mAb	WB,IHC,IF	H,M,R
A0036	STIP1 Rabbit mAb	WB,IHC	H,M,R
A9036	STK3/MST2 Rabbit mAb	WB,IHC,IF	H,M,R
A19803	STK33 Rabbit mAb	WB,IF	H,M
A2275	STK39 Rabbit mAb	WB,IHC,IF	H,M,R
A19243	STX1A Rabbit mAb	WB,IHC,IF	H,M,R
A8334	SULT2A1 Rabbit mAb	WB,IF	H,R
A19121	SUMO1 Rabbit mAb	WB,IHC,IF	H,M,R
A5066	SUMO2/3 Rabbit mAb	WB,IHC,IF	H,M,R
A9016	SUMO4 Rabbit mAb	WB,IHC,IF	H,M,R
A19782	SUN2 Rabbit mAb	WB,IHC	H,M,R
A4838	Survivin Rabbit mAb	WB	M,R
A19542	Synapsin II Rabbit mAb	WB,IF	H,M,R
A19122	Synaptophysin Rabbit mAb	WB,IF	M,R
A3712	Syntaxin 3 Rabbit mAb	WB,IHC	H,M,R
A5996	Syntaxin 4 Rabbit mAb	WB,IF	H,M,R
A19813	Syntaxin Rabbit mAb	WB,IHC,IF	H,M,R
A5497	Syntenin Rabbit mAb	WB,IHC,IF	H,M,R
A3306	TAB1 Rabbit mAb	WB,IF	H,M,R
A19824	TAB3 Rabbit mAb	WB	H,M
A19617	TACC3 Rabbit mAb	WB,IHC,IF	H,M,R

货号	名称	应用	物种
A19810	Talin 2 Rabbit mAb	WB,IHC,IF	H,M,R
A9105	TARBP2 Rabbit mAb	WB,IHC,IF	H,M,R
A3458	TBK1/NAK Rabbit mAb	WB,IF,IP	H,M,R
A19550	TBR1 Rabbit mAb	WB,IF	H,M,R
A19786	TC10/RHOQ Rabbit mAb	WB,IHC	H,M,R
A19548	TCF4/TCF7L2 Rabbit mAb	WB,IHC	H,M,R
A3887	TCP1 alpha Rabbit mAb	WB,IHC	H
A4700	TCP1 beta/CCT2 Rabbit mAb	WB,IHC	H,M,R
A2394	TCTP/TPT1 Rabbit mAb	WB,IHC,IP	H,M,R
A19123	TDP-43/TARDB Rabbit mAb	WB,IHC,IF	H,M,R
A3173	TFAM Rabbit mAb	WB,IHC,IF	H
A0548	TFE3 Rabbit mAb	WB,IHC	H,M,R
A6859	TFF1/pS2 Rabbit mAb	WB,IF	H
A2407	TGF beta induced TGFBI Rabbit mAb	WB,IHC	H,M,R
A19618	TGN46/TGOLN2 Rabbit mAb	WB,IHC,IF	H
A4024	Thioredoxin 1 (Trx1/TXN) Rabbit mAb	WB,IHC	H
A4424	Thioredoxin 2 (Trx2/TXN2) Rabbit mAb	WB,IF	H,M,R
A4298	THOC4/ALYREF Rabbit mAb	WB,IHC	H,M,R
A1536	TIMM50 Rabbit mAb	WB,IHC,IF	H,M,R
A9663	TIRAP Rabbit mAb	WB,IHC	H,M,R
A3528	TLE1 Rabbit mAb	WB,IHC	H,M,R
A19125	TLR2 Rabbit mAb	WB,IF	H
A4116	TMEFF2 Rabbit mAb	WB,IHC	H,R
A19770	TMOD2 Rabbit mAb	WB,IHC,IF	H,M,R
A9126	TMPRSS2 Rabbit mAb	WB,IHC	H,M,R
A19128	TNFAIP3 Rabbit mAb	WB,IHC	H
A19403	TOM20 Rabbit mAb	WB,IHC,IF,IP	H,M,R
A9666	TOMM22 Rabbit mAb	WB,IHC	H,M,R
A4210	tPA/PLAT Rabbit mAb	WB,IF	H,M,R
A4637	TPPP/p25 Rabbit mAb	WB,IHC	M,R
A19130	Transferrin Rabbit mAb	WB,IHC	H
A19774	Transportin 3 (TNPO3) Rabbit mAb	WB,IHC,IF	H,M,R
A20379	TriMethyl-Histone H3-K36 Rabbit mAb	DB,WB,ChIP	H,M,R
A2693	TrkA+B+C Rabbit mAb	WB,IHC	H,M,R

货号	名称	应用	物种
A8723	Tropomyosin 1 Rabbit mAb	WB,IHC,IF,IP	M,R
A4740	Troponin I2 (TNNI2) Rabbit mAb	WB,IF	M,R
A19540	TSC2 Rabbit mAb	WB	H
A3292	TTF1 Rabbit mAb	WB,IHC,IF	H,M,R
AC039	Tubulin alpha-1B/K-alpha-1 Rabbit mAb	WB,IHC,IF	H,M,R
A19805	Tubulin beta-1 chain Rabbit mAb	WB,IHC,IF	H,M,R, Cow,Dog,Monkey
A5955	TWEAKR/Fn14/CD266 Rabbit mAb	WB,IF	H,M,R
A3686	U1A/SNRPA Rabbit mAb	WB,IF	H,M,R
A4552	U2AF2/U2AF65 Rabbit mAb	WB,IHC	H,M,R
A8662	UAP1 Rabbit mAb	WB,IF	H,M,R
A9749	UAP56/DDX39B Rabbit mAb	WB,IHC	H,M,R
A9254	UBA1 Rabbit mAb	WB,IHC,IF	H,M,R
A0503	UBE2B Rabbit mAb	WB,IHC	H,M,R
A5206	UBE2C Rabbit mAb	WB,IF	H
A4396	UBE2I Rabbit mAb	WB,IHC	H,M,R
A9257	Ube2N/Ubc13 Rabbit mAb	WB,IHC,IF	H,M,R
A19686	Ubiquitin Rabbit mAb	WB,IHC,IF	H,M,R
A1372	UCHL3 Rabbit mAb	WB,IHC	H,M,R
A19526	UCP1+UCP3 Rabbit mAb	WB,IHC,IF	M,R
A5959	ULK3 Rabbit mAb	WB,IF	H,M,R
A5071	UPF1/RENT1 Rabbit mAb	WB,IF,IP	H,M,R
A4366	UQCRC2 Rabbit mAb	WB,IHC	H,M,R
A9517	UQCRFS1 Rabbit mAb	WB,IF	H,M,R
A4454	USP10 Rabbit mAb	WB,IF	H,M
A1393	USP39 Rabbit mAb	WB,IHC	H,R
A2389	USP39 Rabbit mAb	WB,IHC,IF	H,M,R
A3448	USP7/HAUSP Rabbit mAb	WB,IHC	H,M,R
A9782	USP9X Rabbit mAb	WB,IHC,IF	H,M,R
A3533	VAMP1 Rabbit mAb	WB,IHC,IF	H,M,R
A4235	VAMP2 Rabbit mAb	WB,IF	H,M,R
A8812	VAMP3/Cellubrevin Rabbit mAb	WB,IF	H,M
A4728	VAMP8 Rabbit mAb	WB,IHC	H,M,R
A3548	VAV2 Rabbit mAb	WB,IHC	H,M

货号	名称	应用	物种
A9583	VAV3 Rabbit mAb	WB,IHC	H,M
A19131	VCAM1 Rabbit mAb	WB,IF	M,R
A1402	VCP Rabbit mAb	WB,IF	H,M,R
A19707	VDAC1 Rabbit mAb	WB,IHC,IF	H,M,R
A19655	Versican Rabbit mAb	WB,IHC,IF	H,M,R
A4691	Villin1 Rabbit mAb	WB,IF	H
A2752	Vinculin Rabbit mAb	WB,IHC	H,M,R
A19798	VPS11 Rabbit mAb	WB,IHC	H,M
A19240	VPS28 Rabbit mAb	WB,IHC	H,M,R
A13523	VWF Rabbit mAb	WB,IHC	H,M
A5923	WAPL Rabbit mAb	WB,IF	H,M
A19601	WAVE2/WASF2 Rabbit mAb	WB,IHC	H,M,R
A19765	WDR4 Rabbit mAb	WB,IHC	H,M,R
A19133	WNT5A Rabbit mAb	WB,IHC	H
A2446	WT1 Rabbit mAb	WB,IHC	M,R
A9022	Xanthine Oxidase (XDH) Rabbit mAb	WB,IHC,IF	H,M,R
A19241	XPD/ERCC2 Rabbit mAb	WB,IHC	H,M,R
A19134	YAP1 Rabbit mAb	WB,IHC,IF	H
A3534	YB-1/YBX1 Rabbit mAb	WB,IHC,IF,IP	H,M,R
A19569	YY1 Rabbit mAb	WB,IF	H,M,R
A9536	ZAP70 Rabbit mAb	WB,IF	H
A19785	ZBTB7A/FBI-1/LRF Rabbit mAb	WB,IHC,IF	H
A11680	Zinc- α 2-glycoprotein (ZAG/AZGP1) Rabbit mAb	WB,IHC	H,R
A1487	ZNF259 Rabbit mAb	WB,IHC	H,M
A2298	Zyxin Rabbit mAb	WB,IF	H,M
A2235	α -Actin-1 (ACTA1) Rabbit mAb	WB,IHC,IF	H,M,R
A2280	α -Actin-1 (ACTA1) Rabbit mAb	WB,IF	H,M,R
A2319	α -Actin-1 (ACTA1) Rabbit mAb	WB,IF	H,M,R
A3379	α -Actinin-4 Rabbit mAb	WB,IF	H,M,R
A17910	α -Smooth Muscle Actin (ACTA2) Rabbit mAb	WB,IHC,IF	H,M,R
A20407	α -Synuclein Rabbit mAb	WB,IHC,IF	H,M,R
A6830	α -Tubulin Rabbit mAb	WB,IHC,IF,IP	H,M,R
A17913	β III-Tubulin Rabbit mAb	WB,IF,IP	H,M,R
A4798	β II-Tubulin/ β 2-Tubulin Rabbit mAb	WB,IHC,IF	H,M,R

货号	名称	应用	物种
A12289	β -Tubulin Rabbit mAb	WB,IF	H,M,R
A4157	γ -Catenin Rabbit mAb	WB,IHC,IF	H,M,R
A9657	γ -Tubulin Rabbit mAb	WB,IHC,IF	H,M,R
A11399	δ -Catenin/p120 Catenin Rabbit mAb	WB,IHC	H,M,R