

HRPT2 / CDC73 Rabbit pAb

Catalog No.: A0636

Basic Information

Observed MW

Refer to figures

Calculated MW

61kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

This gene encodes a tumor suppressor that is involved in transcriptional and post-transcriptional control pathways. The protein is a component of the PAF protein complex, which associates with the RNA polymerase II subunit POLR2A and with a histone methyltransferase complex. This protein appears to facilitate the association of 3' mRNA processing factors with actively-transcribed chromatin. Mutations in this gene have been linked to hyperparathyroidism-jaw tumor syndrome, familial isolated hyperparathyroidism, and parathyroid carcinoma.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

79577

Swiss Prot

Q6P1J9

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 435-505 of human HRPT2 / CDC73 (Q6P1J9).

Synonyms

HYX; FIHP; HPTJT; HRPT1; HRPT2; C1orf28; HRPT2 / CDC73

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from mouse testis, using HRPT2 / CDC73 Rabbit pAb (A0636).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.