

EV71 3C Rabbit pAb

Catalog No.: A10003 **4 Publications**

Basic Information

Observed MW

Refer to Figures

Calculated MW

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Enterovirus

Background

Recommended Dilutions

WB 1:1000 - 1:5000

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Recombinant protein of Enterovirus EV71 3C

Synonyms

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.