

SMPD3 Rabbit pAb

Catalog No.: A10197SP

5 Publications

Basic Information

Observed MW

71 kDa

Calculated MW

71 kDa/70 kDa

Category

Primary antibody

Applications

WB, IP, IHC-P, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to enable phosphatidic acid binding activity; phosphatidylserine binding activity; and sphingomyelin phosphodiesterase activity. Predicted to be involved in positive regulation of exosomal secretion and sphingomyelin metabolic process. Predicted to act upstream of or within several processes, including animal organ development; enzyme linked receptor protein signaling pathway; and sphingolipid metabolic process. Predicted to be located in Golgi apparatus and plasma membrane. Predicted to be active in cytoplasm. Biomarker of pulmonary emphysema.

Recommended Dilutions

WB 1:1000 - 1:5000**IP** 0.5 µg - 4 µg antibody for
400 µg - 600 µg extracts of
whole cells**IHC-P** 1:200 - 1:800**ELISA** Recommended starting
concentration is 1 µg/mL.
Please optimize the
concentration based on your
specific assay
requirements. For high-ratio
antibody dilutions
(≥1:10000) a sequential
dilution method is strongly
recommended to ensure
measurement accuracy.

Immunogen Information

Gene ID

55512

Swiss Prot

Q9NY59

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

NSMASE2

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C


Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Validation Data



Western blot analysis of lysates from SH-SY5Y cells using SMPD3 Rabbit pAb (A10197SP) at 1:2000 dilution incubated overnight at 4°C.


Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25 µg per lane.


Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).


Exposure time: 20 s.



Immunohistochemistry analysis of paraffin-embedded Human small intestine tissue using SMPD3 Rabbit pAb (A10197SP) at a dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Buffer (pH 6.0) prior to IHC staining.



Immunohistochemistry analysis of paraffin-embedded Mouse brain tissue using SMPD3 Rabbit pAb (A10197SP) at a dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Buffer (pH 6.0) prior to IHC staining.



Immunohistochemistry analysis of paraffin-embedded Rat brain tissue using SMPD3 Rabbit pAb (A10197SP) at a dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Buffer (pH 6.0) prior to IHC staining.