

FSTL3 Rabbit pAb

Catalog No.: A10279 **1 Publications**

Basic Information

Observed MW

27kDa

Calculated MW

28kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse, Rat

Background

Follistatin-like 3 is a secreted glycoprotein of the follistatin-module-protein family. It may have a role in leukemogenesis.

Recommended Dilutions

WB 1:1000 - 1:2000

Immunogen Information

Gene ID

10272

Swiss Prot

O95633

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 27-263 of human FSTL3 (NP_005851.1).

Synonyms

FLRG; FSRP; FSTL3

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using FSTL3 antibody (A10279) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 60s.