

CLPTM1L Rabbit pAb

Catalog No.: A10468 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

62kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

The protein encoded by this gene is a membrane protein whose overexpression in cisplatin-sensitive cells causes apoptosis. Polymorphisms in this gene have been reported to increase susceptibility to several cancers, including lung, pancreatic, and breast cancers.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

81037

Swiss Prot

Q96KA5

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 32-280 of human CLPTM1L (NP_110409.2).

Synonyms

CRR9; CLPTM1L

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.