

BABAM1 Rabbit pAb

Catalog No.: A10677

Basic Information

Observed MW

Refer to figures

Calculated MW

37kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Involved in several processes, including mitotic G2 DNA damage checkpoint signaling; positive regulation of DNA repair; and protein K63-linked deubiquitination. Located in cytosol and nuclear body. Part of BRCA1-A complex and BRISC complex.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

29086

Swiss Prot

Q9NWW8

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

NBA1; HSPC142; MERIT40; C19orf62; BABAM1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.