

AHI1 Rabbit pAb

Catalog No.: A10678

Basic Information

Observed MW

Calculated MW

137kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Mouse, Rat

Background

This gene is apparently required for both cerebellar and cortical development in humans. This gene mutations cause specific forms of Joubert syndrome-related disorders. Joubert syndrome (JS) is a recessively inherited developmental brain disorder with several identified causative chromosomal loci. Alternatively spliced transcript variants encoding different isoforms have been identified.

Recommended Dilutions

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

54806

Swiss Prot

Q8N157

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ORF1; AHI-1; JBTS3; dj71N10.1; AHI1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.