

HESX1 Rabbit pAb

Catalog No.: A10696 **1 Publications**

Basic Information

Observed MW

Refer to Figures

Calculated MW

21kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

This gene encodes a conserved homeobox protein that is a transcriptional repressor in the developing forebrain and pituitary gland. Mutations in this gene are associated with septo-optic dysplasia, HESX1-related growth hormone deficiency, and combined pituitary hormone deficiency.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

8820

Swiss Prot

Q9UBX0

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-185 of human HESX1 (NP_003856.1).

Synonyms

ANF; RPX; CPHD5; HESX1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.