

FCRL6 Rabbit pAb

Catalog No.: A10700

Basic Information

Observed MW

Refer to figures

Calculated MW

48kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables MHC class II protein binding activity and protein phosphatase binding activity. Predicted to be involved in cell surface receptor signaling pathway. Located in external side of plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

343413

Swiss Prot

Q6DN72

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

FcRH6; FCRL6

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.