

TMEM106A Rabbit pAb

Catalog No.: A10718

Basic Information

Observed MW

Refer to figures

Calculated MW

29kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to be involved in several processes, including glycoprotein biosynthetic process; positive regulation of cytokine production; and positive regulation of intracellular signal transduction. Predicted to be located in plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

113277

Swiss Prot

Q96A25

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

TMEM106A

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.