

FAM3B Rabbit pAb

Catalog No.: A1082

1 Publications

Basic Information

Observed MW

Calculated MW

26kDa

Category

Primary antibody

Applications

ELISA, IHC-P, IF/ICC

Cross-Reactivity

Human, Mouse, Rat

Background

Involved in insulin secretion. Located in extracellular exosome.

Recommended Dilutions

IHC-P 1:100 - 1:200

IF/ICC 1:100 - 1:200

Immunogen Information

Gene ID

54097

Swiss Prot

P58499

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 30-235 of human FAM3B (NP_478066.3).

Synonyms

2-21; ORF9; PANDER; PRED44; C21orf11; C21orf76; FAM3B

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.


Validation Data




Immunofluorescence analysis of paraffin-embedded mouse pancreas using FAM3B Rabbit pAb (A1082) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunohistochemistry analysis of paraffin-embedded human colon using FAM3B Rabbit pAb (A1082) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma using FAM3B Rabbit pAb (A1082) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry analysis of paraffin-embedded human stomach using FAM3B Rabbit pAb (A1082) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.