

# PIN1b Rabbit pAb

Catalog No.: A10989

## Basic Information

---

### Observed MW

Refer to Figures

### Calculated MW

64kDa

### Category

Primary antibody

### Applications

WB, ELISA

### Cross-Reactivity

Oryza sativa

## Background

---

## Recommended Dilutions

---

**WB** 1:500 - 1:2000

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

---

### Gene ID

4330700

### Swiss Prot

### Immunogen

Recombinant protein of Oryza sativa LOC4330700

### Synonyms

PIN1; REH1; OsPIN1; OsPIN1a; OsPIN1b; OsJ\_08350; OJ1734\_E02.14; PIN1b

## Contact

---

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.