

ENHO Rabbit pAb

Catalog No.: A11733

Basic Information

Observed MW

Calculated MW

8kDa

Category

Primary antibody

Applications

IHC-P, ELISA

Cross-Reactivity

Mouse, Rat

Background

Predicted to enable hormone activity. Predicted to be involved in positive regulation of Notch signaling pathway. Predicted to be located in extracellular region. Predicted to be active in plasma membrane.

Recommended Dilutions

IHC-P 1:50 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

375704

Swiss Prot

Q6UWT2

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

UNQ470; C9orf165; ENHO

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG


Purification

Affinity purification


Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.



Immunohistochemistry analysis of paraffin-embedded Rat testis using ENHO Rabbit pAb (A11733) at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to IHC staining.



Immunohistochemistry analysis of paraffin-embedded Mouse testis using ENHO Rabbit pAb (A11733) at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to IHC staining.