

Tyro3 Rabbit pAb

Catalog No.: A12241 **2 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

97kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

The gene is part of a 3-member transmembrane receptor kinase receptor family with a processed pseudogene distal on chromosome 15. The encoded protein is activated by the products of the growth arrest-specific gene 6 and protein S genes and is involved in controlling cell survival and proliferation, spermatogenesis, immunoregulation and phagocytosis. The encoded protein has also been identified as a cell entry factor for Ebola and Marburg viruses.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

7301

Swiss Prot

Q06418

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 791-890 of human Tyro3 (NP_006284.2).

Synonyms

BYK; Dtk; RSE; Rek; Sky; Tif; Etk-2; Tyro3

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.