

# SCGB2A1 Rabbit pAb

Catalog No.: A12406 **1 Publications**

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

11kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human

## Background

Involved in androgen receptor signaling pathway. Located in extracellular space.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

4246

### Swiss Prot

O75556

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human SCGB2A1 (NP\_002398.1).

### Synonyms

LPHC; LPNC; MGB2; UGB3; SCGB2A1

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.