

NDP Rabbit pAb

Catalog No.: A12470

Basic Information

Observed MW

Calculated MW

15kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Rat

Background

This gene encodes a secreted protein with a cystein-knot motif that activates the Wnt/beta-catenin pathway. The protein forms disulfide-linked oligomers in the extracellular matrix. Mutations in this gene result in Norrie disease and X-linked exudative vitreoretinopathy.

Recommended Dilutions

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

4693

Swiss Prot

Q00604

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ND; EVR2; FEVR; NDP

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.