

# RAD51C Rabbit pAb

Catalog No.: A12487

## Basic Information

### Observed MW

### Calculated MW

42kDa

### Category

Primary antibody

### Applications

ELISA,IF/ICC

### Cross-Reactivity

Human

## Background

This gene is a member of the RAD51 family. RAD51 family members are highly similar to bacterial RecA and *Saccharomyces cerevisiae* Rad51 and are known to be involved in the homologous recombination and repair of DNA. This protein can interact with other RAD51 paralogs and is reported to be important for Holliday junction resolution. Mutations in this gene are associated with Fanconi anemia-like syndrome. This gene is one of four localized to a region of chromosome 17q23 where amplification occurs frequently in breast tumors. Overexpression of the four genes during amplification has been observed and suggests a possible role in tumor progression. Alternative splicing results in multiple transcript variants.

## Recommended Dilutions

IF/ICC 1:50 - 1:100

## Immunogen Information

### Gene ID

5889

### Swiss Prot

O43502

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human RAD51C (NP\_002867.1).

### Synonyms

FANCO; R51H3; BROVCA3; RAD51L2; RAD51C

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data

---



Immunofluorescence analysis of MCF7 cells using RAD51C antibody (A12487).