

# ABI Rabbit pAb

Catalog No.: A12786 **3 Publications**

## Basic Information

### Observed MW

Refer to Figures

### Calculated MW

73kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Zea mays

## Background

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

542523

### Swiss Prot

### Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of zea mays VP1 (NP\_001105540.1).

### Synonyms

vp1; GRMZM2G133398; ABI

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.