

# NRARP Rabbit pAb

Catalog No.: A12875

## Basic Information

### Observed MW

12kDa

### Calculated MW

12kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Mouse

## Background

Involved in Notch signaling pathway. Acts upstream of or within negative regulation of Notch signaling pathway and positive regulation of canonical Wnt signaling pathway.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

441478

### Swiss Prot

Q7Z6K4

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-114 of human NRARP (NP\_001004354.1).

### Synonyms

NRARP

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

## Validation Data

---



Western blot analysis of lysates from mouse kidney, using NRARP Rabbit pAb (A12875) at 1:3000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25 $\mu$ g per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Enhanced Kit (RM00021).  
Exposure time: 90s.