

# THAP3 Rabbit pAb

Catalog No.: A12899

## Basic Information

### Observed MW

27kDa

### Calculated MW

27kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Mouse, Rat

## Background

Predicted to enable DNA binding activity and metal ion binding activity. Involved in positive regulation of transcription by RNA polymerase II.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

90326

### Swiss Prot

Q8WTV1

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-239 of human THAP3 (NP\_001182682.1).

### Synonyms

THAP3

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

## Validation Data

---



Western blot analysis of various lysates using THAP3 Rabbit pAb (A12899) at 1:3000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 90s.