

AVEN Rabbit pAb

Catalog No.: A12910 **1 Publications**

Basic Information

Observed MW

55kDa

Calculated MW

39kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Rat

Background

Involved in negative regulation of apoptotic process. Located in membrane.

Recommended Dilutions

WB 1:1000 - 1:3000

Immunogen Information

Gene ID

57099

Swiss Prot

Q9NQ51

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 203-362 of human AVEN (NP_065104.1).

Synonyms

PDCD12; AVEN

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using AVEN Rabbit pAb (A12910) at 1:3000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.