

# MRGPRX1 Rabbit pAb

Catalog No.: A13142

## Basic Information

### Observed MW

36kDa

### Calculated MW

36kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Mouse, Rat

## Background

Enables transmembrane signaling receptor activity. Involved in cell surface receptor signaling pathway and response to chloroquine. Predicted to be located in cell surface.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

259249

### Swiss Prot

Q96LB2

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 263-322 of human MRGPRX1 (NP\_671732.3).

### Synonyms

GPCR; MGRG2; MRGX1; SNSR4; MRGPRX1

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

## Validation Data

---



Western blot analysis of various lysates using MRGPRX1 Rabbit pAb (A13142) at 1:3000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 90s.