

CCAR1 Rabbit pAb

Catalog No.: A13595

Basic Information

Observed MW

132kDa

Calculated MW

133kDa

Category

Primary antibody

Applications

ELISA, WB, IP

Cross-Reactivity

Human, Mouse, Rat

Background

Enables RNA polymerase II cis-regulatory region sequence-specific DNA binding activity; nuclear receptor coactivator activity; and transcription corepressor activity. Involved in positive regulation of cell migration and positive regulation of cell population proliferation. Acts upstream of or within positive regulation of apoptotic process. Located in nuclear envelope lumen.

Recommended Dilutions

WB 1:500 - 1:1000

IP 0.5µg-4µg antibody for
200µg-400µg extracts of
whole cells

Immunogen Information

Gene ID

55749

Swiss Prot

Q8IX12

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CCAR1 (NP_060707.2).

Synonyms

CCAR1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Immunoprecipitation analysis of 200 µg extracts of HeLa cells using 3 µg CCAR1 antibody (A13595). Western blot was performed from the immunoprecipitate using CCAR1 antibody (A13595) at a dilution of 1:500.



Western blot analysis of various lysates using CCAR1 Rabbit pAb (A13595) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 5s.