

# RFFL Rabbit pAb

Catalog No.: A13604

## Basic Information

### Observed MW

### Calculated MW

41kDa

### Category

Primary antibody

### Applications

ELISA

### Cross-Reactivity

Human

## Background

Enables enzyme binding activity; p53 binding activity; and ubiquitin protein ligase activity. Involved in cellular protein metabolic process; negative regulation of cysteine-type endopeptidase activity involved in execution phase of apoptosis; and negative regulation of signal transduction. Located in endosome membrane and plasma membrane.

## Recommended Dilutions

### ELISA

Recommended starting concentration is 1  $\mu\text{g/mL}$ . Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

117584

### Swiss Prot

Q8WZ73

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

CARP2; FRING; CARP-2; RNF189; RNF34L; RIFIFYLIN; RFFL

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at  $-20^{\circ}\text{C}$ . Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.