

ZNF396 Rabbit pAb

Catalog No.: A13836

Basic Information

Observed MW

41kDa

Calculated MW

39kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Enables protein heterodimerization activity; protein homodimerization activity; and sequence-specific double-stranded DNA binding activity. Involved in negative regulation of transcription, DNA-templated. Located in cytoplasm and nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

252884

Swiss Prot

Q96N95

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ZSCAN14; ZNF396

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using ZNF396 Rabbit pAb (A13836) at 1:3000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 50s.