

ROBO4 Rabbit pAb

Catalog No.: A14134

Basic Information

Observed MW

Refer to figures

Calculated MW

107kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Predicted to enable cell-cell adhesion mediator activity. Involved in angiogenesis and establishment of endothelial barrier. Located in extracellular exosome. Implicated in aortic valve disease 3.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

54538

Swiss Prot

Q8WZ75

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

MRB; AOVD3; ECSM4; ROBO4

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.