

PLCXD2 Rabbit pAb

Catalog No.: A14164

Basic Information

Observed MW

Calculated MW
35kDa

Category
Primary antibody

Applications
WB,IF/ICC,ELISA

Cross-Reactivity
Human, Mouse, Rat

Background

Predicted to enable phosphoric diester hydrolase activity. Predicted to be involved in lipid catabolic process and signal transduction. Predicted to be located in nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

IF/ICC 1:50 - 1:100

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID 257068	Swiss Prot Q0VAA5
--------------------------	-----------------------------

Immunogen
Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms
PLCXD2


Contact

		400-999-6126
		cn.market@abclonal.com.cn
		www.abclonal.com.cn

Product Information

Source Rabbit	Isotype IgG	Purification Affinity purification
Storage Store at -20°C		

Validation Data



Immunofluorescence analysis of A549 cells using PLCXD2 Rabbit pAb (A14164). Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution.