

RPUSD2 Rabbit pAb

Catalog No.: A14397

Basic Information

Observed MW

Refer to Figures

Calculated MW

61kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse, Rat

Background

Enables pseudouridine synthase activity. Involved in mRNA pseudouridine synthesis.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

27079

Swiss Prot

Q8IZ73

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

PUS9; C18B11; C15orf19; RPUSD2

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.