

CDK17 Rabbit pAb

Catalog No.: A14761

Basic Information

Observed MW

Calculated MW

60kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Mouse, Rat

Background

The protein encoded by this gene belongs to the cdc2/cdkx subfamily of the ser/thr family of protein kinases. It has similarity to a rat protein that is thought to play a role in terminally differentiated neurons. Alternatively spliced transcript variants encoding different isoforms have been found.

Recommended Dilutions

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

5128

Swiss Prot

Q00537

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

PCTK2; PCTAIRE2; CDK17

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.