

ZNF273 Rabbit pAb

Catalog No.: A14846

Basic Information

Observed MW

Calculated MW

65kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Rat

Background

This gene is a member of the krueppel C2H2-type zinc-finger protein family and encodes a protein with 13 C2H2-type zinc fingers and a KRAB domain. This nuclear protein is involved in transcriptional regulation. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

10793

Swiss Prot

Q14593

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

HZF9; ZNF273

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.