

VPS26B Rabbit pAb

Catalog No.: A15205

Basic Information

Observed MW

Calculated MW
39kDa

Category
Primary antibody

Applications
ELISA

Cross-Reactivity
Human

Background

Predicted to be involved in intracellular protein transport and retrograde transport, endosome to Golgi. Predicted to act upstream of or within cellular response to interferon-gamma. Predicted to be located in early endosome and late endosome. Predicted to be part of retromer complex. Predicted to be active in endosome.

Recommended Dilutions

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID	Swiss Prot
112936	Q4G0F5

Immunogen
Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms
Pep8b; VPS26B

Contact

	400-999-6126
	cn.market@abclonal.com.cn
	www.abclonal.com.cn

Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage
Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.