

# THNSL1 Rabbit pAb

Catalog No.: A15515

## Basic Information

### Observed MW

83kDa

### Calculated MW

83kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human, Mouse, Rat

## Background

THNSL1 (Threonine Synthase Like 1) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include pyridoxal phosphate binding and threonine synthase activity. An important paralog of this gene is THNSL2.

## Recommended Dilutions

WB 1:200 - 1:2000

## Immunogen Information

### Gene ID

79896

### Swiss Prot

Q8IYQ7

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 534-743 of human THNSL1 (NP\_079114.3).

### Synonyms

TSH1; THNSL1

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

## Validation Data

---



Western blot analysis of various lysates using THNSL1 Rabbit pAb (A15515) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 30s.