

DNER Rabbit pAb

Catalog No.: A15551

Basic Information

Observed MW

80kDa

Calculated MW

78kDa

Category

Primary antibody

Applications

WB, IF/ICC, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to enable Notch binding activity. Involved in central nervous system development. Located in dendrite; early endosome; and plasma membrane.

Recommended Dilutions

WB 1:200 - 1:2000

IF/ICC 1:50 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

92737

Swiss Prot

Q8NFT8

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

bet; UNQ26; DNER

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C

Validation Data

