

DDX51 Rabbit pAb

Catalog No.: A15580

Basic Information

Observed MW

72kDa

Calculated MW

72kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse, Rat

Background

Enables RNA binding activity. Predicted to be involved in rRNA processing. Located in membrane.

Recommended Dilutions

WB 1:200 - 1:2000

Immunogen Information

Gene ID

317781

Swiss Prot

Q8N8A6

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 367-666 of human DDX51 (NP_778236.2).

Synonyms

DDX51

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using DDX51 Rabbit pAb (A15580) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 5s.