

# PRR11 Rabbit pAb

Catalog No.: A15599

## Basic Information

---

### Observed MW

45kDa

### Calculated MW

40kDa

### Category

Primary antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human, Mouse, Rat

## Background

---

Involved in regulation of cell cycle. Located in cytoplasm and nucleus.

## Recommended Dilutions

---

**WB** 1:500 - 1:2000

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

---

### Gene ID

55771

### Swiss Prot

Q96HE9

### Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

### Synonyms

PRR11

## Contact

---

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.