

MPHOSPH6 Rabbit pAb

Catalog No.: A15775

Basic Information

Observed MW

22kDa,23kDa

Calculated MW

19kDa

Category

Primary antibody

Applications

ELISA,WB

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to enable RNA binding activity. Involved in maturation of 5.8S rRNA. Located in cytosol; nucleolus; and nucleoplasm. Colocalizes with nuclear exosome (RNase complex).

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

10200

Swiss Prot

Q99547

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human MPHOSPH6 (NP_005783.2).

Synonyms

MPP; MPP6; MPP-6; MPHOSPH6

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates using MPHOSPH6 Rabbit pAb (A15775) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 10s.