

MED31 Rabbit pAb

Catalog No.: A15831

Basic Information

Observed MW

16kDa

Calculated MW

16kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to enable transcription coregulator activity and ubiquitin protein ligase activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to act upstream of or within limb development and negative regulation of fibroblast proliferation. Predicted to be located in nucleoplasm. Predicted to be part of core mediator complex and mediator complex.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

51003

Swiss Prot

Q9Y3C7

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Soh1; CGI-125; 3110004H13Rik; MED31

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using MED31 Rabbit pAb (A15831) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 5s.