

SH2D1B Rabbit pAb

Catalog No.: A16010

Basic Information

Observed MW

Calculated MW

15kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Human

Background

By binding phosphotyrosines through its free SRC (MIM 190090) homology-2 (SH2) domain, EAT2 regulates signal transduction through receptors expressed on the surface of antigen-presenting cells (Morra et al., 2001 [PubMed 11689425]).

Recommended Dilutions

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

117157

Swiss Prot

O14796

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

EAT2; SH2D1B

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.