

PHF2 Rabbit pAb

Catalog No.: A16063

Basic Information

Observed MW

Refer to figures

Calculated MW

121kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Human, Rat

Background

This gene encodes a protein which contains a zinc finger-like PHD (plant homeodomain) finger, distinct from other classes of zinc finger motifs, and a hydrophobic and highly conserved domain. The PHD finger shows the typical Cys4-His-Cys3 arrangement. PHD finger genes are thought to belong to a diverse group of transcriptional regulators possibly affecting eukaryotic gene expression by influencing chromatin structure.

Recommended Dilutions

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

5253

Swiss Prot

O75151

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

GRC5; KDM7C; JHDM1E; CENP-35; PHF2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.